# RIDASCREEN®

## Chlamydia trachomatis IgA (K 2911)

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Čeština

# RIDASCREEN® Chlamydia trachomatis IgA (K 2911)

For *in vitro* diagnostic use. This test is an enzyme immunoassay for the determination of IgA antibodies against Chlamydia trachomatis in human serum.

## Test procedure

1. Bring the microwell plate [Plate] and reagents to room temperature (20-25°C).
2. Dilute the wash buffer [Wash] with distilled water 1:20.
3. Place enough wells in the frame plate [Plate] for Positive control [Control IgA +], Negative control [Control IgA -], Cut-off control [Cut Off IgA] and samples. The cut-off control must be carried out in duplicate.
4. Dilute the controls and the serum samples in the plate [Plate] 1:21 with the sample buffer [Diluent]. Absorb the samples in the [Plate] (e.g. with RIDA® RF-Absorbens, Article no. Z 0202) and do not adjust to the required dilution until this has been done.  
First place 50 µl RIDA® RF-Absorbens in the wells for the samples  
Add 5 µl serum and mix  
Place 5 µl controls in the corresponding wells  
Add 100 µl [Diluent] to the controls  
Add 50 µl [Diluent] to the samples and mix
5. Incubate at 37°C for 45 minutes.
6. Empty the microwell plate and then wash it 5 times with 300 µl diluted wash buffer.
7. Place 100 µl conjugate [Conjugate IgA] in each well.
8. Incubate at 37°C for 30 minutes.
9. Empty the microwell plate [Plate] and then wash it 5 times with 300 µl diluted wash buffer.
10. Place 100 µl [Substrate] in each well
11. Incubate at room temperature (20-25 °C) for 20 minutes.
12. Place 50 µl stop reagent [Stop] in each well
13. Carry out a photometric measurement at 450/620 nm.

The test has been carried out correctly when the following conditions have been fulfilled

	OD
Negative control	< 0.55
Positive control	> 0.9
Cut-off control (average value)	< 0.7 x OD for the positive control > 1.5 x OD for the negative control

## Evaluation

1. The average absorbance is calculated for the cut-off control.
2. The sample index is obtained by dividing the absorbance for the sample by the calculated average value.

For example:

Cut-off control 1	OD = 0.821
Cut-off control 2	OD = 0.865
Average value	= 0.843
Sample	OD = 1.508

$$\text{Sample index} = \frac{1.508}{0.843} = 1.79$$

Evaluating the sample index




negative	equivocal	positive
< 0.9	0.9-1.1	> 1.1

## Explanation of symbols

### General symbols:

	For <i>in vitro</i> diagnostic use
	Consult instructions* for use
	Lot number
	Expiry
	Store at
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

### Test-specific symbols:

	Microwell plate
	Sample dilution buffer
	Wash buffer
	Positive control IgA
	Negative control IgA
	Cut-off control IgA
	Anti-human IgA conjugate
	Substrate
	Stop reagent

\* Please obtain instructions from [www.r-biopharm.com](http://www.r-biopharm.com) or enquire from your local R-Biopharm agent.

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**R-Biopharm AG**