

**Calibrex™** *digital*  
your universal preference

# bottle-top dispensers



Designed for performance handling of liquids from a large variety of bottles and flasks, the Calibrex™ instruments combine the latest in dosing technology, high-tech materials and ergonomic design. As a result, users benefit from ease of operation, superior safety and low maintenance.

**The Calibrex™ line features:**

- Superior chemical resistance
- Long-lasting performance stability
- Comfortable and convenient to use
- Solid, yet simple construction
- Instant volume setting
- In-lab calibration
- Autoclavable at 121°C fully assembled

**The Calibrex™ models**  
**520 521**

 **SOCOREX**  
SWISS

# Calibrex™ digital

your universal preference



## ► Friendly digital volume setting

Adjustment is fast and precise even with a gloved hand. An efficient click-stop mechanism prevents unwanted alteration. Large display is easy to read; the window adjusts to desired side of instrument body.

A true digital system contributes to the long-term calibration stability. Each step on the cylindrical cam is pre-calibrated and corresponds to one division on the volume adjustment scale.




## ► Optimised working position

The instrument rotates 360° for adequate positioning on top of the bottle.



## ► Air filter

Air inlet can be enlarged to receive a membrane filter, if airborne contamination is a concern.



## ► Extension tubing

Autoclavable PTFE delivery Jet-Pen and extension tubing help dispense into vessels with maximum comfort within a 60 cm distance. This accessory fits Calibrex™ dispensers up to 20 mL.

## Calibrex™ 520 up to 10 mL

Three colour-coded models with adjustable volumes are ranging 0.25 to 2 mL, 1 to 5 mL and 1 to 10 mL. Fixed versions also available in any volume. The Calibrex™ combine a clever, yet simple construction and a universal chemical compatibility. Among other features, they are particularly appreciated for their smooth, effortless activation and their sturdiness.



520

0.25 - 2 mL  
1 - 5 mL  
1 - 10 mL

## Absolute stability

Stand holds dispenser up to 50 mL. Most useful with small-size bottles.



## Bulk container aspiration

Work station facilitates liquid intake from drum or other remote container (<10 m distance, <2 m elevation).





## Calibrex™ 521 up to 100 mL

Three 521 sizes expand the Calibrex™ line with larger dispensing volumes ranging from 2 to 20 mL, 10 to 50 mL and 10 to 100 mL. Similar to the 520 model in many respects such as construction quality, sturdiness, chemical resistance, smooth plunger travel, in-lab calibration, autoclaving, the 521 range has a distinct black design. Also fitted with the Solid Calibration adjustment system, the Calibrex™ 521 has an additional fine-tuning ring.

## 521



2 - 20 mL  
10 - 50 mL  
10 - 100 mL

### Performance

Volume mL	Inaccuracy (E %)*		Imprecision (CV%)	
	Min. vol.	Max. vol.	Min. vol.	Max. vol.
<b>Calibrex™ 520</b>				
0.25 - 2	< ± 3.0 %	< ± 0.6 %	< 0.5 %	< 0.1 %
1 - 5	< ± 2.0 %	< ± 0.6 %	< 0.5 %	< 0.1 %
1 - 10	< ± 1.5 %	< ± 0.6 %	< 0.5 %	< 0.1 %
<b>Calibrex™ 521</b>				
2 - 20	< ± 1.5 %	< ± 0.6 %	< 0.5 %	< 0.1 %
10 - 50	< ± 1.5 %	< ± 0.6 %	< 0.5 %	< 0.1 %
10 - 100	< ± 1.5 %	< ± 0.6 %	< 0.5 %	< 0.1 %

\*Fixed volume models: performance may vary according to desired volume



### ► In-lab calibration

The Calibrex™ dispensers are factory calibrated. Both 520 and 521 models can be easily recalibrated by means of a dedicated adjustment screw.



### ► Easy maintenance

Disassembling/reassembling is facilitated by the limited number of well-connected parts - no tool needed. Key elements such as valves, plunger, barrel and delivery jet can be removed in seconds for cleaning. All Calibrex™ dispensers are autoclavable fully assembled at 121°C.



### ► Superior material selection

Special attention is paid to component materials, providing for long instrument life. Parts coming in contact with the liquid flow are chemically inert.

Parts	520 model	521 model
Feed tube	PTFE	
Intake valve	Ceramic	Borosilicate glass
Valve balls	Pyrex glass	Ceramic, ruby
Valve spring	Platinum-iridium	
Barrel	Borosilicate glass	
Barrel plate	PTFE	
Plunger	Glass with PFA coating	
Outlet valve	Ceramic	
Body	ETFE	
Delivery jet assembly	PTFE/ETFE	



### ► Bottle-neck adapters

Both 520 and 521 models have a 32 mm base thread. They come with three bottle-neck adapters in 28, 40 and 45 mm diameters (520) or 38, 40, 45 mm diameters (521), bringing to four the number of suited bottle diameters. See p.4 for bottles and accessories.

# bottle-top dispensers

## Ordering information – Instruments

Standard supply includes instrument, 3 adapters (see right column on p.3), feed tube 300 mm (375 mm with 521.100), QC certificate and operating instructions.

Volume mL	Increment mL	Cat. No. Adj. volume	Cat. No. Fixed volume*	Volume mL	Increment mL	Cat. No. Adj. volume	Cat. No. Fixed volume*
<b>Calibrex™ 520</b>				<b>Calibrex™ 521</b>			
0.25 – 2	0.05	520.002	520.F02	2 – 20	0.1	521.020	521.F020
1 – 5	0.1	520.005	520.F05	10 – 50	1.0	521.050	521.F050
1 – 10	0.25	520.010	520.F10	10 – 100	1.0	521.100	521.F100

\*Specify desired fixed volume within instrument range when ordering

## Ordering information – Accessories

### Work station

Description	Cat. No.
Stand for dispenser up to 50 mL	320.SB050
Stand for remote liquid intake, dispenser up to 50 mL*	320.BC050
Stand for remote liquid intake, dispenser 100 mL*	320.BC100

\*Feed tube not included, to be ordered separately

### PTFE tubings

Description	Length	Int. Ø	Cat. No.
Extention tubing for dispenser up to 20 mL	600 mm	2.2 mm	1.524
Replacement feed tubing, 520 all sizes	300 mm	5 mm	511.707
Replacement feed tubing 521, 20 mL	300 mm	7.5 mm	521.621
Replacement feed tubing 521, 50 mL	300 mm	10 mm	521.622
Replacement feed tubing 521, 100 mL	375 mm	10 mm	521.623
Feed tubing by the meter, uncut	specify	5 mm	511.709
Feed tubing by the meter, uncut	specify	7.5 mm	521.705
Feed tubing by the meter, uncut	specify	10 mm	521.706

### Adapters – Replacement and additional adapters of various kind

Type	Ext. Ø of bottle neck	Cat. No.
Threaded, PP	22, 25, 28, 30, 34, 36, 38, 40, 45 mm	GLP + Ø
Threaded, PTFE	36, 38, 40, 45 mm	GLT + Ø
Tapered, PP	18.8, 24, 29.2, 45 mm	NSP + Ø
Tapered, PTFE	18.8, 24, 29.2 mm	NST + Ø

### Bottles – Large selection of bottles supplied with screw cap

Material	Volume mL	Neck, ext. Ø	Cat. No.	Material	Volume mL	Neck, ext. Ø	Cat. No.
Amber glass	100	32 mm	314.0100	Clear borosilicate glass	250	45 mm	314.0250C
	250	32 mm	314.0250		500	45 mm	314.0500C
	500	32 mm	314.0500		1000	45 mm	314.1000C
	1000	45 mm	314.1000		2000	45 mm	314.2000C
	2500	45 mm	314.2500		500	32 mm	314.0500P
Amber glass, with handle	2500	45 mm	314.2500H	Clear Pyrex glass, with connection neck	250	25 mm	315.0250
PE coated amber glass	500	32 mm	314.0500PE		500	25 mm	315.0500
	1000	45 mm	314.1000PE		1000	32 mm	315.1000
	2500	45 mm	314.2500PE		2500	45 mm	315.2500

## Your local distributor

## Safety first

A protection sleeve prevents barrel breakage and eliminates any risk of spilling. Cap-on delivery jet prevents liquid loss. Both models have been approved for safety (GS tested) by the German TÜV organisation.



## QC and warranty

Each Calibrex™ dispenser bears its own serial number and passes a strict quality control attested by an individual certificate.

Performance values obtained by a smooth and constant plunger movement, with Aqua bidest. at constant temperature ( $\pm 0.5^\circ\text{C}$ ) comprised between 20 and 25°C, according to EN ISO 8655/DIN 12650.

One-year warranty



**SOCOREX ISBA SA**  
**Champ-Colomb 7, P.O. Box**  
**1024 Ecublens/Lausanne**  
**Switzerland**  
**socorex@socorex.com**  
**www.socorex.com**  
**Tel: +41 (0)21 651 6000**  
**Fax: +41 (0)21 651 6001**