## **RIDASCREEN<sup>®</sup> Diphtherie IgG**

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Enzyme immunoassay for the detection of  
IgG antibodies against diphtheria toxoid in serum

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RIDA® and RIDASCREEN®  
are registered trademarks of R-Biopharm GmbH  
Manufacturer: R-Biopharm GmbH, Darmstadt, Germany

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**Contents**

	page
1. Intended use.....	14
2. General.....	14
3. Test principle.....	15
4. Reagents provided.....	16
5. Materials required but not provided.....	16
6. Warnings and precautions for the users.....	17
7. Storage instructions.....	17
8. Indication of instability or deterioration of reagents.....	18
9. Specimen collection and storage.....	18
10. Test procedure.....	18
11. Analysis.....	21
12. Remarks about the test procedure and interpretation.....	23

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## RIDASCREEN® Diphtherie IgG

Enzyme immunoassay for the detection of IgG antibodies against diphtheria toxoid in serum

### 1. Intended use

The RIDASCREEN® Diphtherie test is an enzyme immunoassay (EIA) for the detection of IgG antibodies against diphtheria toxoid in human serum.

### 2. General

The symptoms of diphtheria are caused by *Corynebacterium diphtheriae*. This is a gram-positive, rod-shaped bacteria which is often thicker at the ends. The length is 2 - 6 µm and the diameter is 0.5 - 1 µm. Typically, the rods are arranged in V, Y or palisade formation. A characteristic feature during selective culture on tellurite-containing media is the blackening of the colonies caused by the reduction of the telluric salt to metallic tellurium.

The pathogenic behavior of *C. diphtheriae* is caused by an exotoxin, the gene for which is located in a prophage. The toxin is a 62 kilodalton polypeptide composed of two fragments of different functionality. The toxin bonds to receptors on the target cell via Fragment B and is absorbed into the cell. Once inside the cell, Fragment A blocks the synthesis of protein irreversibly.

Two courses of the infection can be distinguished for diphtheria: local infection and primary systemic intoxication. In the case of local infection, the mucous membranes of the mouth and throat are predominantly infected. After an incubation time of two to five days, uncharacteristic symptoms develop initially, such as headache, exhaustion, difficulty in swallowing and moderate fever. The throat reddens and flat pseudomembranes develop as a result of the lesions at the site of infection. These membranes may spread to the larynx resulting in a croup, with differing degrees of shortness of breath, and choking attacks in the most severe cases.

It is rarely found that the nasal mucous membranes, conjunctivas or vulva are infected during the course of local infection. Navel diphtheria can develop in the case of new-born children and wound diphtheria in the case of adults. Systemic intoxication results in fulminant diphtheria. The local symptoms in the throat are significantly more pronounced. Due to the effect of the toxin on peripheral nerves and heart, symptoms of paralysis develop and, in the most severe cases, circulatory distress. Cardiac death may occur suddenly as a result of damage to the cardiac muscle (early death on the fifth or sixth day of the disease). Heart failure can also be expected during the recovery period.

Therapy involves administration of an antiserum against the toxin by injection as early as possible and antibiotics. Administering antibiotics on their own, however, is not sufficient. When treating the local infection by administering the antiserum by injection, improvement is rapid in most cases. However, for systemic intoxication fatalities are high in spite of administering the antiserum.

The diagnosis is established from clinical symptoms and via cultural evidence of the pathogen.

Diphtheria is a rare disease owing to the consistent vaccination of infants. Vaccination is carried out using a toxoid obtained by treating the toxin with formaldehyde. This toxoid has the same antigenic properties as the toxin without the toxic effects. IgG vaccine titers, even the immunoglobulins present after recovery from an infection, can therefore be found in a toxoid-based ELISA, such as the RIDASCREEN® Diphtheria IgG EIA. Comparability of the measured values is ensured by matching to international units.

### 3. Test principle

Diphtheria toxoid is bound on the surface of the microtiter wells. Diluted serum samples and controls are pipetted into the wells and incubated at 37 °C in a humid chamber. Present antibodies bind to the immobilized antigen. Unbound material is removed in a washing step.

In a second step, a peroxidase-conjugated anti-human antibody (anti-IgG) is added. After incubation, unbound conjugate is removed by washing. Substrate (H<sub>2</sub>O<sub>2</sub>/TMB) is added to the wells and incubated at 37 °C in a humid chamber. The enzyme bound in the wells converts the colorless substrate to a blue color. Addition of stop solution converts the color from blue to yellow. The absorption is measured at 450 nm wavelength (reference wavelength ≥ 620 nm).

#### 4. Reagents provided

The reagents in one package are sufficient for 96 determinations.  
Each test kit contains:

- 1x 12 Microwell Strips (breakable) with 8 wells each in a frame;  
coated with diphtheria toxoid;  
in a resealable foil bag
- 1 x RIDASCREEN® SeroPP (110 ml);  
Sample Diluent, ready to use;  
contains 0.01 % Thimerosal
- 1 x RIDASCREEN® SeroWP (100 ml, 10x conc.);  
Washing Buffer, contains 0.2 % Bronidox-L
- 1 x Standard Control (2.5 ml);  
diluted human serum, ready to use;  
contains 0.01 % Thimerosal
- 1 x Negative Control (1.2 ml);  
diluted human serum, ready to use;  
contains 0.01 % Thimerosal
- 1 x RIDASCREEN® SeroG HD (12 ml);  
Anti-human-IgG-Conjugate;  
peroxidase-conjugated antibody (rabbit);  
ready to use
- 1 x RIDASCREEN® SeroSC (12 ml);  
H<sub>2</sub>O<sub>2</sub>/tetramethylbenzidine (TMB), ready to use
- 1 x RIDASCREEN® SeroStopp (12 ml);  
1 M sulfuric acid
- 1 x Instructions for use

#### 5. Reagents required but not provided

##### 5.1 Reagents

- Distilled or deionized water

##### 5.2 Accessories

- Test tubes
- Vortex mixer
- Micropipets for volumes of 10 - 100 µl and 1 ml
- Microplate washer or multichannel pipet
- Microplate reader (450 nm, reference wavelength ≥ 620 nm)
- Absorbent paper

#### 6. Warnings and precautions for the users

The control sera (standard control and negative control) have been tested for HIV- and HCV-Ab as well as for HBsAg and were found to be negative. However, they as well as the patient samples should be considered potentially contagious and be treated with the necessary safety precautions.

The standard control and the negative control as well as the sample diluent contain 0.01 % Thimerosal. Contact with skin or mucous membranes must be avoided.

The washing buffer contains Bronidox-L. Contact with skin or mucous membranes must be avoided.

Hydrogen peroxide can cause cauterization. Handle with care!

The stop solution contains 1 N sulfuric acid. Avoid contact with skin and clothing!

All reagents and materials coming in contact with potential infectious specimens must be treated with disinfectants or autoclaved at 121 °C for at least one hour.

An exchange of the controls between kits of different lot numbers is not possible.

#### 7. Storage instructions

All reagents have to be stored at 2 - 8 °C and can be used up to the expiry date printed on the labels. Microbial contamination has to be avoided. A quality warranty cannot be given beyond the kit expiration date.

The diluted washing buffer (RIDASCREEN® SeroWP) can be used up to the expiry date printed on the labels if stored at 2 - 8 °C.

Allow reagents and microwell strips to get room temperature before use. To avoid moisture within the strips, do not take the strips out of the foil bag before having reached room temperature. The foil bag should be opened with a pair of scissors without detaching the fastener. Return any unused strips to the foil bag, reseal and store them directly at 2 - 8 °C.

It is absolutely necessary to avoid contamination of the substrate (RIDASCREEN® SeroSC) with the conjugate solution, since this results in a coloration of the substrate. In the same way, the colorless substrate must be protected from exposure to direct light to avoid deterioration or coloration by autoxidation. If the substrate has turned blue, the reagent should be discarded.

## 8. Indication of instability or deterioration of reagents

The following criteria may indicate a reagent deterioration:

- a turbidity or a blue coloration of the colorless substrate prior to its use
- an absorbance value of the negative control at 450/620 nm > 0.3
- an absorbance value of the positive control at 450/620 nm out of the values given in the attached tables of values

## 9. Specimen collection and storage

The RIDASCREEN® Diphtherie IgG EIA has been evaluated for the investigation of human serum samples. Repeated freezing and thawing of the samples as well as microbial contamination must be avoided. The application of heat treated, lipemic, hemolytic, icteric or turbid samples can lead to wrong results.

The sample material can be stored for up to 1 week at 2 - 8 °C if the test cannot be carried out immediately. A prolonged storage of samples is possible at -20 °C.

Diluted samples can be used for up to 7 hours if they have been stored at 2 - 8 °C.

## 10. Test procedure

### 10.1. Preliminary comments

The RIDASCREEN® EIAs for serological antibody detection contain some test independent reagents. These reagents are labeled RIDASCREEN® Sero (SeroWP, SeroPP, SeroG HD, SeroSC, SeroStopp). They are not kit specific and can be changed between single test kits with correspondingly labeled reagents. But, this is not valid for the controls.

Bring all reagents and the microwell strips to room temperature before use. Mix the reagents well before use. Reproducibility in any EIA depends on exact pipetting, the observance of incubation times and temperature and the consistency of wash sequences. A deviation from the given incubation times and temperatures leads to an excessive difference of the standard control serum to the given target value, which is no longer covered by the determined area of values. During the washing steps, take care that all wells are filled with buffer and that the liquid is completely removed from the wells. Do not allow microwells to dry between steps.

Avoid direct sunlight during all incubations. Covering the microtiter plate is recommended.

Except the washing buffer, all reagents are ready to use.

### 10.2. Preparation of the washing buffer

1 part of the concentrated washing buffer (RIDASCREEN® SeroWP) is diluted with 9 parts of distilled water. Crystals in the buffer concentrate can be dissolved in a waterbath at 37 °C. Add 100 ml of the concentrated washing buffer to a

1000 ml graduated cylinder. Bring the final volume to 1000 ml with distilled or deionized water. The diluted washing buffer can be used up to the expiry date printed on the labels if stored at 2 - 8 °C.

### 10.3. Preparation of the samples

Before starting the test, serum samples have to be diluted with the sample diluent (RIDASCREEN® SeroPP).

Dilution of the serum sample 1:100:

e.g. 10 µl serum + 990 µl sample diluent

#### **Attention!**

**Negative control and standard control are ready to use and must not be diluted.**

### 10.4. First incubation

After insertion of a sufficient number of cavities into the microwell holder, 100 µl of the diluted sera and ready to use controls are each added to the corresponding wells, the position A1 remains empty (blank). One negative control and two standard controls are carried along. The plate is incubated in a humid chamber at 37 °C for 30 min. The bottom of the cavities should not be in touch with materials that conduct temperature well (metal or wet paper).

A1	blank
B1	negative control
C1	standard control
D1	standard control
E1, F1	patient serum 1, 2, etc.

#### **Attention!**

**The ELISA plate should not be placed in a cold incubation container which is heated to 37 °C during the incubation. The container must already be adapted to 37 °C prior to the incubation.**

### 10.5. Washing

Decant or aspirate all microwells into a waste container with a disinfectant. Ensure complete removal of the liquid from the microwells by tapping the inverted plate onto absorbent paper. Then wash all wells 4 times with 300 µl of prepared washing buffer. Be sure to remove residual washing solution by firmly tapping the inverted microwells on absorbent paper after single washing steps.

**If a microplate washer is used, take care that the washer is adjusted to the used microplate type.**

## 10.6. Second incubation

Add 100 µl of RIDASCREEN® SeroG HD (anti-*human-IgG*-conjugate) into each well (including A1). Incubate the plate for 30 min at 37 °C in a humid chamber.

## 10.7. Washing

Wash 4 times according to step 10.5.

## 10.8. Third incubation

Add 100 µl of RIDASCREEN® SeroSC into each well. Incubate the plate for 30 min at 37 °C in a humid chamber. Following the incubation, the reaction is stopped by adding 100 µl of RIDASCREEN® SeroStopp to each well. After careful mixing (soft tapping on the edge of the plate) the absorbance is measured in a microplate reader at 450 nm (reference wavelength ≥ 620 nm). Correct for blank value against position A1.

### Remark:

**To remove moisture, wipe the bottom of the microplate before measuring.**

### Summary of the test procedure

1. Bring all reagents to room temperature
2. Dilute the RIDASCREEN® SeroWP
3. Dilute the serum samples
4. Pipet 100 µl of the standard control, the negative control or the diluted samples into the microwells; 30 minutes incubation at 37 °C in a humid chamber
5. Discard the incubate and wash 4 times with 300 µl of washing buffer
6. Add 100 µl of RIDASCREEN® SeroG HD; 30 minutes incubation at 37 °C in a humid chamber
7. Discard the incubate and wash 4 times with 300 µl of washing buffer
8. Add 100 µl of RIDASCREEN® SeroSC; 30 minutes incubation at 37 °C in a humid chamber
9. After addition of 100 µl RIDASCREEN® SeroStopp spectrophotometric determination at 450 nm (reference wavelength ≥ 620 nm)

## 11. Analysis

Analysis of the test can be done in three different ways:

1. using the standard curve enclosed
2. using the table of values (see data sheet enclosed)
3. mathematically by the 4-parameter calculation or the  $\alpha$ -method

**The blank (A1) must be subtracted from all OD values prior to the analysis.**

### 11.1. Quality control

For the quality control, standard control (double determination) and negative control must be carried along with each test procedure. The test was carried out correctly if the mean absorbance value of the standard control at 450/620 nm corresponds to the value range of the data sheet enclosed. If the single measurements diverge from the mean absorbance value more than 20 %, the test must be repeated. The negative control must show an absorbance value at 450/620 nm < 0.3.

### 11.2. Analysis with the standard curve

If the analysis is to be carried out with the standard curve, first the daily variation must be corrected by a double determination of the standard control (see also pt. 11.1.). The correction factor F is calculated with the given target value of the standard control and the actually measured standard control value. The Lot-dependent target value is given on the data sheet enclosed.

$$F = \frac{\text{target value of the standard control}}{\text{mean absorbance value of the standard control}}$$

All OD values of the samples are multiplied by the factor F. With these corrected values, the corresponding IUnit value is read from the standard curve.

### 11.3. Analysis with the table of values

The mean absorbance value of the standard control (see also pt. 11.1.) determines the column of the value range in the attached value table, which is decisive for the actual measurement. Within this column the measured absorbance value of the sample is classed to the suited absorbance range and the corresponding titer in IUnits/ml is read in the second column from the left.

IUnits/ml	range of standard control	
	0.77 - 0.81	
< 0.10	< 0.10	
0.10 - 0.29	0.10 - 0.22	
0.30 - 0.49	0.23 - 0.35	
0.50 - 0.99	0.36 - 0.64	
1.00 - 1.49	0.65 - 0.91	
1.50 - 2.00	0.92 - 1.13	
> 2.00	> 1.13	

Fig. 1: Example (Excerpt from a Lot dependent data sheet)

The mean absorbance value of the standard control in a measurement is 0.81 for example. In this case, the column with the value range from 0.77 to 0.81 is decisive for result determination. Then, a patient sample with an absorbance of 0.41 lies in a titer range from 0.50 to 0.99 IUnits/ml. (The foregoing values are an example and can differ from the actual values of the Lot-dependent data sheet.)

#### 11.4. Mathematical analysis

Needed values for a mathematical analysis by the 4-parameter calculation or the  $\alpha$ -method are given on the data sheet enclosed.

#### 11.5. Test result

IU/ml	recommended immunization
< 0.1	initial immunization
0.1 - 0.9	booster immunization
1.0 - 1.4	booster immunization in 5 years
1.5 - 2.0	booster immunization in 7 years
> 2.0	booster immunization in 10 years

## 12. Remarks about the test procedure and interpretation

The information under Section 11.5. is oriented towards a publication by Pietsch in 1993. The author assumed that immunization could start with a titer of 0.1 IU/ml (international units per ml), but would only be sufficient with titers greater than 1.0 IU/ml. Basic immunization should be carried out in the case of titers smaller than 0.1 IU/ml.

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