**INTENDED USE**

The CerTest FOB-Transferrin is a one step coloured chromatographic immunoassay for the qualitative determination of human hemoglobin (Hb) and human transferrin (Tf) in stool samples to detect gastrointestinal bleeding.

**INTRODUCTION**

Colorectal cancer is a leading cause of illness and death in the Western world. Screening with fecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps, will bleed, for which may be detected by an occult blood test. When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact hemoglobin, intact heme, and heme-derived porphyrins in amounts that depend on the site, the amount of bleeding and the transit time through the gut. Immunochemical tests detect intact or nearly intact human hemoglobin, being a very specific technique to detect the loss of blood from the lower intestine, because blood from lower sites is less degraded during transit. The hemoglobin is unstable in feces, leading to false negative results. The detection of fecal transferrin, which is more stable than hemoglobin, provides an alternative way of diagnosing the disease in the upper digestive tract. Transferrin is a blood-derived component that may be leaked into gastrointestinal tract and then discharged with the feces in gastrointestinal bleeding diseases. Transferrin is stable in feces and a good marker to detect loss of blood from the upper and lower intestine (gastrointestinal bleeding). This immunochromatographic assay detects human hemoglobin and human transferrin in stool samples simultaneously, by obtaining more accurate testing results.

**PRINCIPLE OF THE TEST**

The CerTest FOB-Transferrin is a qualitative immunochromatographic assay for the determination of human hemoglobin and human transferrin in stool samples. The membrane is pre-coated with mouse monoclonal antibodies on both test bands (result region), against human hemoglobin and human transferrin.

During testing, the sample is allowed to react with the coloured conjugates (anti-human hemoglobin antibodies-blue microspheres and anti-human transferrin antibodies-red microspheres) pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. Different coloured bands will be visible, depending upon the hemoglobin or the transferrin content of the sample. These bands are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a green coloured band always appears. The presence of this green band serves as 1) verification that sufficient volume is added, 2) proper flow obtained and 3) an internal control for the reagents.

**STORAGE AND STABILITY**

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

**PRECAUTIONS**


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered as potentially hazardous and handled as if they were infectious agents.
- The tests should be discarded in a proper biohazard container after testing.

**SPECIMEN COLLECTION AND PREPARATION**

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

**Specimen preparation (see illustration):**

- (1) Unscrew the tap and use the stick by introducing four times into the fecal specimen to pick up a little sample (approx. 150mg), if the stool sample were liquid take 150 µL by using a pipette and adding the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

**Precautions:** Patients suffering from menstrual period, bleeding haemorrhoids, blood in urine or strain during bowel movement should not collect samples.

**MATERIALS PROVIDED**

- Device tests
- Instructions for use
- Stool collection tubes-sample diluent


**MATERIALS REQUIRED BUT NO PROVIDED**

- Specimen collection container
- Disposable gloves
- Timer


**TEST PROCEDURE**

**Allow the test, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.**

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the top (3).
2. Remove the CerTest FOB-Transferrin card device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100 µL into the circular window marked with the arrow (4).
4. Read the result at **10 minutes** (the coloured bands appear). Do not read the test result after more than 10 minutes.

**INTERPRETATION OF RESULTS (PLEASE REFER TO THE ILLUSTRATION BELOW)**

**NEGATIVE:** Only one GREEN band appears across the central window in the site marked with the letter C (control line). No occult blood presents: neither human haemoglobin nor human transferrin.

**Tf and Hb POSITIVE:** In addition to the GREEN control band, a RED band (Tf test line) and a BLUE one (Hb test line) also appear in the site marked with the letter T (result region). That probably would mean a lower gastrointestinal bleeding disease (colorectal cancer).

**Tf POSITIVE:** In addition to the GREEN control band, a RED band (Tf test line) also appears in the site marked with the letter T (result region). That probably would mean an upper gastrointestinal bleeding disease (the human haemoglobin was probably degraded in the gastrointestinal tract).

**Hb POSITIVE:** In addition to the GREEN control band, a BLUE band (Hb test line) also appears in the site marked with the letter T (result region). That probably would mean a lower gastrointestinal bleeding disease (the concentration of human transferrin in blood is 100 times less than human haemoglobin, therefore, Hb positive means not much blood present in feces).

**INVALID:** A total absence of the control coloured band (GREEN) regardless the appearance or not of the result lines (RED or BLUE). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are likely the reasons for control line failure. Review the procedure and repeat the tests using a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

**NOTES ON THE INTERPRETATION OF RESULTS**

The intensity of the RED/BLUE coloured band in the result region (T) will vary depending on the concentration of human hemoglobin or human transferrin in the specimen. However, neither the quantitative value, nor the rate of increase in hemoglobin or transferrin can be determined by this qualitative test.

**QUALITY CONTROL**

Internal procedural controls are included in the test. A GREEN line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

**LIMITATIONS**

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of stool sample could result in wrong results (brown bands appear or absence of the control coloured band).
3. Patients suffering from menstrual period, bleeding hemorrhoids, blood in urine or strain during bowel movement should not collect samples.
4. Positive results confirm the presence of human hemoglobin or/and human transferrin in fecal samples; nevertheless, it can be also due to several causes besides colorectal bleeding, such as hemorrhoids, blood in urine or stomach irritations. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood in the stool.
5. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.

**PERFORMANCE**

**Sensitivity**

A sample containing human hemoglobin at concentration equal to or higher than 30 ng/mL and human transferrin at concentration equal to or higher than 4 ng/mL produces positive results when using CerTest FOB-Transferrin test.

Different hemoglobin and transferrin dilutions were tested directly in the extraction buffer or spiked in a negative stool sample in accordance with the kit instructions to determine the detection limit of the test.

The detection of human hemoglobin with CerTest FOB-Transferrin test showed >99% of sensitivity compared to a commercial guaiac assay.

The detection of human hemoglobin in stool samples with CerTest FOB-Transferrin test showed >99 % of sensitivity compared with others commercial rapid tests (ImmunoTech OccuTech (YD Diagnostics) and Human Hexagon OBTI).

**Specificity**

The CerTest FOB-Transferrin test is specific for human hemoglobin and human transferrin, showing no cross-reaction with hemoglobin or transferrin from bovine and pig.

The detection of human haemoglobin with CerTest FOB-Transferrin Test showed >99% of specificity compared to a commercial guaiac assay, >99% compared with other commercial rapid tests (ImmunoTech OccuTech (YD Diagnostics) and Human Hexagon OBTI).

The use of mouse monoclonal antibodies in the elaboration of CerTest FOB-Transferrin test assures high degree of specificity to detect human haemoglobin and transferrin.

**REFERENCES**

1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.
4. Chien-Hua Chiang et al., "A Comparative Study of Three Fecal Occult Blood Test in Upper Gastrointestinal Bleeding". Kaohsiung J. Med. Sci. May 2006, 22; 223-228.
5. Virtanen et al., "Higher concentrations of serum transferrin receptor in children than in adults", Am. J. Clin. Nutr. 1999; 69: 256-60.

**SYMBOLS FOR IVD COMPONENTS AND REAGENTS**

	<i>In vitro</i> diagnostic device		Lot code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Do not use if package damaged

**CerTest**  
BIOTEC S.L.  
María de Luna 11  
E-50018 Zaragoza (SPAIN)  
www.certest.es

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1. Towler BP, Irwig J, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Faecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.
4. Chien-Hua Chiang et al., "A Comparative Study of Three Faecal Occult Blood Test in Upper Gastrointestinal Bleeding". Kaohsiung J. Med. Sci. May 2006, 22; 223-228.
5. Virtanen et al., "Higher concentrations of serum transferrin receptor in children than in adults", Am. J. Clin. Nutr. 1999; 69: 256-60.

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