INTENDED USE

The CerTest HAV is a one step coloured chromatographic immunoassay for the qualitative detection of Hepatitis A virus in stool samples.

INTRODUCTION

Hepatitis A is a liver infection caused by the hepatitis A virus (HAV). The virus is spread when an uninfected (or unvaccinated) person eats or drinks something contaminated by the stool of an HAV-infected person: this is called faecal-oral transmission. The symptoms of hepatitis A range from mild to severe, and can include fever, malaise, loss of appetite, diarrhoea, nausea, abdominal discomfort, dark-colored urine and jaundice (a yellowing of the skin and whites of the eyes). Not everyone who is infected will have all of the symptoms. Adults have signs and symptoms of illness more often than children, and the severity of disease and mortality increases in older age groups. Infected children under six years of age do not usually experience noticeable symptoms, and only 10% develop jaundice. Among older children and adults, infection usually causes more severe symptoms, with jaundice occurring in more than 70% of cases. Most people recover in several weeks - or sometimes months - without complications.

PRINCIPLE OF THE TEST

The CerTest HAV is a qualitative immunochromatographic assay for the determination of Hepatitis A virus in fecal samples. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, against viral antigens. During testing, the sample is allowed to react with the coloured conjugate (anti-Hepatitis A virus mouse monoclonal antibodies-red microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugate. The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers and the assay should be done right after collection. The samples can be stored in the refrigerator (2-4°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):

- (1) Unscrew the cap and use the stick to pick up a little sample, if the stool sample was liquid take 100 µL using a pipette, and add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

MATERIALS PROVIDED

- Device tests (CerTest HAV)
- Instructions for use
- Stool collection tubes with sample diluent


MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


TEST PROCEDURE

Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until ready to perform the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion. Cut the end of the cap (3).
2. Remove the CerTest HAV device from its sealed bag just before using.
3. Use a separate stool collection tube and device for each sample or control. Dispense 4 drops or 100 µL into the circular window marked with an arrow (4).
- 4.- Read the result at **15 minutes** (the coloured bands appear).

INTERPRETATION OF RESULTS (please refer to the illustration below)

NEGATIVE: Only one GREEN band appears across the central window in the site marked with the letter C (control line). A negative result indicates no HAV infection. A confirmed diagnosis should be made after all clinical and laboratory findings.

POSITIVE: In addition to the GREEN control band, a RED band (test line) also appears in the site marked with the letter T (result line). This result indicates a recent and contagious infection with HAV.

INVALID: A total absence of the control coloured band (GREEN) regardless the appearance or not of the result line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. A green line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear).
3. After one month of infection, the number of viruses in feces is decreasing, making the sample less reactive. Stool samples should be collected previously to the onset of symptoms to stop the spread of viruses.
4. This test provides a presumptive diagnosis for Hepatitis A virus infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

EXPECTED VALUES

Negative results are expected in healthy infants and young children, also in healthy adults.

PERFORMANCE

Sensitivity and specificity

It was studied some stool samples from patients with the same hepatitis A symptoms: fever, diarrhoea and dark-colored urine and some patients with jaundice. All samples were confirmed as positive results by HAV-Antigen EIA assay (Mediagnost®). It was performed an evaluation using CerTest HAV and the results showed >99% of sensitivity and >99% of specificity. The antibodies used to elaborate this test recognise epitopes present in antigens found in stool patients.

REFERENCES

1. COULEPIS, A. et al. "Detection of Hepatitis A Virus and Antibody by Solid-Phase Radioimmunoassay and Enzyme-Linked Immunosorbent Assay with monoclonal antibodies". *Journal of Clinical Microbiology*, July 1985, Vol. 22, No. 1, p. 119-124.
2. BOROVEC, S. et al. "Synthesis and Assembly of Hepatitis A Virus-Specific Proteins in BS-C-1 Cells". *Journal of Virology*. June 1993, Vol. 67, No.6, p. 3095-3102.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Lot code
	Consult instructions for use		Catalogue number
	Keep dry		Contains sufficient for <n> tests
	Temperature limitation		Manufacturer
	Use by		Sample diluent

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