VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

H. pylori + Clarithromycin resistance

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ENGLISH

1. Intended use

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit is designed for the specific identification of *H. pylori* and detection of clarithromycin (CLR) resistance in gastric tissue biopsies from patients with signs and symptoms of gastrointestinal infection. This test is intended to be used as an aid in the diagnosis of *H. pylori* and resistance to clarithromycin in combination with clinical and epidemiological risk factors. DNA is extracted from gastric tissue biopsies, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *H. pylori* and clarithromycin. The assay for detecting CLR resistance is based on detection of point mutations in the 23S rRNA gene.

2. Summary and Explanation

The genus *Helicobacter* belongs to the family *Helicobacteriaceae* of the order *Campylobacterales*. *Helicobacter pylori* (*H. pylori*) is a gram-negative microaerophilic spiral-shaped bacterium, which is able to colonize the mucus layer of the human stomach and the upper part of small intestine (duodenum).

More than half of the world population is estimated to be infected with *H. pylori*, but most individuals are asymptomatic. The new infections might be due to direct person-to-person transmission, via either an oral-oral, fecal-oral route or both. Nevertheless, the role of coccoid form of *H. pylori* as a vehicle of infection from food and water sources should not be discarded. *H. pylori* is involved in the pathogenesis of atrophic gastritis, gastroduodenal ulcer, gastric cancer and gastric mucosa-associated lymphoid tissue (MALT) lymphoma.

At present, several diagnostic assays for *H. pylori* detection are available and grouped as "invasive" or "noninvasive", but none of them can be considered as gold standard alone. Invasive methods include histology, culture and rapid urease testing, which require gastric biopsy specimens obtained by gastroduodenoscopy. Noninvasive approaches include fecal antigen detection, serologic testing, and urea breath testing among others.

Urease is an important factor for the maintenance and virulence of the bacterium in the gastric mucosa. It is composed of two structural subunits encoded by genes, *ureA* and *ureB*, which have been frequently employed as target genes for the specific detection of *Helicobacter pylori*.

Clarithromycin is a bacteriostatic antibiotic mostly used in childhood to treat upper and lower respiratory tract infections, but, its application to treat *H. pylori* is the most used indication. The main action mode of clarithromycin as one of the wide-spectrum antibiotics used in *H. pylori* therapy is to prevent protein translation. Following the first exposure to the clarithromycin, spontaneously mutations (in both 23S rRNA operons) confer *H. pylori* resistance genotype and phenotype. The direct impact of these mutations is emergence of *H. pylori* strains resistant to clarithromycin. To now, two major mutations A2142G and A2143G were listed as main cause of antibiotic resistance in clinical isolates.

3. Principle of the procedure

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit is designed for the diagnosis of *H. pylori*, Clarithromycin resistance and Clarithromycin wild-type sequence in the 23S rRNA in gastric tissue biopsies. After DNA isolation, the identification of *H. pylori*, Clarithromycin resistance and/or Clarithromycin wild-type sequence in the 23S rRNA is performed by the amplification of a conserved region of the *ureB* and 23S rRNA genes respectively, using specific primers and a fluorescent-labeled probe.

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. Point mutations in the 23S rRNA gene of *H. pylori* (A2142G and A2143G), which confer resistance to Clarithromycin are amplified and detected in FAM channel, *H. pylori* DNA targets are amplified and detected in ROX channel, Clarithromycin wild-type sequence in the 23S rRNA gene DNA targets are amplified and detected in HEX, VIC or JOE channel and the internal control (IC) in Cy5 channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>H. pylori</i> + Clarithromycin resistance 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>H. pylori</i> + Clarithromycin resistance Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit with Ref. VS-CLA106L, VS-CLA106H, VS-CLA112L and VS-CLA112H.

Reagent/Material	Description	Color	Amount
<i>H. pylori</i> + Clarithromycin resistance 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>H. pylori</i> + Clarithromycin resistance Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit with Ref VS-CLA113L and VS-CLA113H.

Reagent/Material	Description	Colour	Amount
<i>H. pylori</i> + Clarithromycin resistance 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>H. pylori</i> + Clarithromycin resistance Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit with Ref. VS-CLA172 and VS-CLA172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
 - Vortex.
 - Micropipettes (0.5-20 µL, 20-200 µL).
 - Filter tips.
 - Powder-free disposable gloves.
 - Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-CLA113L, VS-CLA113H, VS-CLA136 and VS-CLA172). Remove any excess air in the pouches prior to sealing.

- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-CLA136 and VS-CLA172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Sample preparation

Gastric tissue biopsies should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples must be frozen at -20°C or -80°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise gastric tissue biopsies as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

8.1.1. DNA extraction

For DNA extraction from gastric tissue biopsies you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended*.
- Invisorb® Spin Universal Kit (Strattec).
- Maxwell®RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).

- ZP02004 MagPurix Tissue DNA Extraction Kit and ZP02006 MagPurix Bacterial DNA Extraction Kit (biopsies specimens), using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.).

*In order to improve the yield and quality of bacterial DNA from biopsies specimens, we recommend the sample pretreatment with Lysozyme at 37°C as described on "Instructions for use" for bacterial DNA isolation. In addition, usage of small elution volumes may raise the DNA/RNA concentration.

8.2. Lyophilized positive control

H. pylori + Clarithromycin resistance Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *H. pylori* + Clarithromycin resistance Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *H. pylori* + Clarithromycin resistance Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	63°C

Table 5. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (Clarithromycin resistance), ROX (*H. pylori*), HEX, JOE or VIC (Clarithromycin wild-type sequence in the 23S rRNA) and Cy5 channels (IC). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *H. pylori* + Clarithromycin resistance in the positive control well. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer’s instructions. To consider that a sample contains *H. pylori* resistant and/or sensitive to clarithromycin it must be *H. pylori* positive. In addition, it is recommended to check the fluorescence level in FAM and HEX channels in order to discriminate between Clarithromycin resistant and/or wild-type sequences.

Using the following table read and analyze the results:

Clarithromycin resistance (FAM)	<i>H. pylori</i> (ROX)	Clarithromycin wild-type sequence (HEX)	Internal control (Cy5)	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>H. pylori</i> positive, Clarithromycin resistance positive and Clarithromycin wild-type sequences in the 23S rRNA positive ^a
+	+	-	+/-	-	+	<i>H. pylori</i> positive and Clarithromycin resistance positive; Clarithromycin wild-type sequence in the 23S rRNA negative ^b
-	+	+	+/-	-	+	<i>H. pylori</i> positive and Clarithromycin wild-type sequence in the 23S rRNA positive, Clarithromycin resistance negative ^c
+	-	+	+/-	-	+	<i>H. pylori</i> negative; Clarithromycin resistance positive and Clarithromycin wild-type sequences in the 23S rRNA positive
-	-	-	+	-	+	<i>H. pylori</i> negative, Clarithromycin resistance negative and Clarithromycin wild-type sequences in the 23S rRNA negative
-	-	+	+/-	-	+	Clarithromycin wild-type sequence in the 23S rRNA positive
+	-	-	+/-	-	+	Clarithromycin resistance positive
-	+	-	+/-	-	+	Experiment fail
-	-	-	-	-	-	Experiment fail
+	+	+	+	+	+	Experiment fail

Table 6. Sample interpretation
 +: Amplification curve
 -: No amplification curve

^a Presence of Clarithromycin resistant strains and wild-type simultaneously in the same clinical sample. Similar fluorescent levels in HEX and FAM channels should be seen. On the contrary, if differences are noted between the fluorescent signal from both channels look at the possibility of interpretation ^b or ^c.

^b Small level of fluorescence can be observed in HEX channels.

^c Small level of fluorescence can be observed in FAM channels.

Figure 1 and 2. Correct run of positive *H. pylori* and clarithromycin resistance sample and small level of fluorescence in HEX channels (wild-type) run on the Bio-Rad CFX96™ Real-Time PCR System (Figure 1) and Correct run of positive *H. pylori* and wild-type sample and small level of fluorescence in FAM channels (clarithromycin resistance) run on the Bio-Rad CFX96™ Real-Time PCR System (Figure 2).


Figure 1



Figure 2


A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 3. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR System.

Negative control

Positive control

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from gastric tissue biopsies.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *H. pylori*, Clarithromycin resistance and Clarithromycin wild-type sequences in the 23S rRNA, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Mutations or polymorphisms in primer or probe binding regions may affect detection of new variants resulting in a false negative result with the VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit.
- Other species different from *H. pylori* could carry different alleles in the 23S rRNA gene and give positive results to the Clarithromycin wild-type sequences and/or the Clarithromycin resistance.

11. Quality control

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit was tested using human gastric tissue biopsies from symptomatic patients. These results were compared with those obtained with a molecular detection method (RIDA®GENE *Helicobacter pylori* (R-Biopharm)).

The results were as follows:

VIASURE <i>H. pylori</i> + Clarithromycin resistance Real Time PCR Detection Kit	RIDA®GENE <i>Helicobacter pylori</i> (R-Biopharm)			
		+	-	Total
	+	99	3*	102
	-	2*	102	104
	Total	101	105	206

Table 7. Comparative results for *H. pylori*.

*The low amount of template DNA in these samples are below the detection limit of the method used.

VIASURE <i>H. pylori</i> + Clarithromycin resistance Real Time PCR Detection Kit	RIDA®GENE <i>Helicobacter pylori</i> (R-Biopharm)			
		+	-	Total
	+	23	0	23
	-	2#	77	79
	Total	25	77	102

Table 8. Comparative results for Clarithromycin resistance.

#These 2 samples are sensitive strains to Clarithromycin according to VIASURE and resistant strains to Clarithromycin according to R-biopharm. They are sensitive strains to Clarithromycin according to sequencing, confirming our results.

VIASURE <i>H. pylori</i> + Clarithromycin resistance Real Time PCR Detection Kit	RIDA®GENE <i>Helicobacter pylori</i> (R-Biopharm)			
		+	-	Total
	+	77	2#	79
	-	0	23	23
	Total	77	25	102

Table 9. Comparative results for Clarithromycin wild-type sequences in the 23S rRNA.

#These 2 samples are sensitive strains to Clarithromycin according to VIASURE and resistant strains to Clarithromycin according to R-biopharm. They are sensitive strains to Clarithromycin according to sequencing, confirming our results.

In conclusion, the results show a high sensitivity and specificity to detect Clarithromycin resistance and *H. pylori* using VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction for Clarithromycin resistance and *H. pylori* (Figure 2, 3 and 4).

Figure 4. Dilution series of Clarithromycin resistance (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR System (channel FAM).


Figure 5. Dilution series of *H.pylori* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR System (channel ROX).



Figure 6. Dilution series of Clarithromycin wild-type sequence (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR System (channel HEX).

12.3. Analytical specificity

The specificity of the *H. pylori* + Clarithromycin resistance assay was confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine. No cross-reactivity was detected against any of the following microorganisms tested, with the exception of *H. heilmannii*, which together with *H. pylori* can cause chronic gastritis in humans.

Cross-reactivity testing					
<i>Helicobacter felis</i> *	-	<i>Campylobacter upsaliensis</i>	-	<i>Candida albicans</i>	-
<i>Helicobacter hepaticus</i> *	-	<i>Campylobacter hyointestinalis</i> *	-	<i>Arcobacter butzleri</i>	-
<i>Helicobacter cinaedi</i> *	-	<i>Proteus vulgaris</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Helicobacter heilmannii</i> *	-/+	<i>Citrobacter freundii</i>	-	<i>Enterococcus faecalis</i>	-
<i>Shigella flexneri</i>	-	<i>Staphylococcus aureus subsp. aureus</i>	-	<i>Yersinia enterocolitica</i> O:3*	-
<i>Shigella dysenteriae</i>	-	<i>Serratia liquefaciens</i>	-	<i>Yersinia enterocolitica</i> O:9	-
<i>Salmonella typhi</i>	-	<i>Vibrio parahaemolyticus</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella paratyphi</i> A	-	<i>Clostridium difficile</i>	-	Adenovirus serotypes 40/41	-
<i>Salmonella paratyphi</i> B	-	<i>Clostridium difficile</i> 027*	-	Adenovirus serotypes 1-5	-
<i>Salmonella typhimurium</i>	-	<i>Clostridium perfringens</i>	-	Adenovirus serotype 8/15/31	-
<i>Salmonella bongori</i>	-	Enterotoxigenic <i>Escherichia coli</i>	-	Rotavirus A	-
<i>Salmonella enteritidis</i>	-	Enterohemorrhagic <i>Escherichia coli</i>	-	Norovirus Genotypes I and II	-
<i>Salmonella enterica subsp. enterica</i>	-	Enteroinvasive <i>Escherichia coli</i>	-	Astrovirus Genotype I-VIII	-
<i>Salmonella pullorum</i>	-	Enteropathogenic <i>Escherichia coli</i>	-	Sapovirus	-
<i>Salmonella gallinarum</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Entamoeba histolytica</i>	-
<i>Campylobacter lari</i> *	-	<i>Listeria monocytogenes</i>	-	<i>Entamoeba dispar</i>	-
<i>Campylobacter fetus</i> *	-	<i>Aeromonas caviae</i>	-	<i>Cryptosporidium parvum/hominis</i>	-
<i>Campylobacter coli</i> *	-	<i>Aeromonas hydrophila subsp. Hydrophila</i> *	-	<i>Giardia intestinalis</i> *	-
<i>Campylobacter jejuni subsp. jejuni</i>	-	<i>Blastocystis hominis</i>	-	<i>Dientamoeba fragilis</i>	-

Table 10. Reference pathogenic microorganisms used in this study.

*These microorganisms showed amplification signal in Clarithromycin wild-type sequences in the 23S rRNA channel.

The specificity of the clarithromycin resistance assay was confirmed by testing a panel consisting of different antimicrobial resistant organism. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing			
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) N315	-	<i>Citrobacter braakii</i> with VIM-1 gene	-
Methicillin-resistant <i>Staphylococcus aureus</i> ST398	-	<i>Citrobacter freundii</i> -complex isolate with KPC-3 and VIM-4 genes	-
mecC Methicillin-resistant <i>Staphylococcus aureus</i>	-	<i>H. pylori</i> Clarithromycin resistant (23S rRNA A2146G)	-/+
cMRSA isolate (oxa ^R , PVL-positive, spa:t 310)	-	<i>H. pylori</i> Clarithromycin resistant (23S rRNA A2147G)	-/+
<i>Serratia marcescens</i> isolate with OXA-48 gene	-	VanA-type <i>Enterococcus avium</i>	-
<i>Klebsiella pneumonia</i> isolate with SHV-1 (non-ESBL), KPC-3, and OXA-48 genes	-	VanA- type <i>Enterococcus faecium</i> LMG16165 strain	-
<i>Klebsiella pneumonia</i> isolate with TEM-1 (non-ESBL), SHV-1 (non-ESBL), CTX-M-2 (ESBL), and KPC-2 genes	-	VanA- type <i>Enterococcus faecium</i> IOWA 1	-
<i>Escherichia coli</i> with OXA-244 gene	-	VanB- type <i>Enterococcus faecium</i> IOWA 2	-
<i>Escherichia coli</i> with TEM-1 (non-ESBL) and IMP-1 genes	-	VanA-type <i>Enterococcus faecalis</i>	-
<i>Enterobacter cloacae</i> with SHV-12 (ESBL), CTX-M-9 (ESBL) and OXA-48 genes	-	VanB-type <i>E. faecalis</i> (Andrewes and Horder) Schleifer and Kilpper-Balz	-
<i>Enterobacter cloacae</i> with TEM-1 (non-ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL) and NDM-1 genes	-	VanC and VanB- types <i>Enterococcus gallinarum</i> ENT20120142	-
<i>Enterobacter cloacae</i> -complex with a NDM-7 gene	-	VanC type- <i>Enterococcus gallinarum</i> (Bridge and Sneath 1982) Collins, Jones, Farrow, Kilpper-Balz and Schleifer 1984 VP	-

Table 11. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit for *Helicobacter pylori* was evaluated against *Helicobacter pylori* strain J99 and Sydney, showing positive result.

The reactivity of VIASURE *H. pylori* + Clarithromycin resistance Real Time PCR Detection Kit for Clarithromycin resistance was evaluated against *H. pylori* strains wild-type and harbour point mutations A2142G or A2143G, showing negative and positives results, respectively.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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

13. Bibliography/Bibliografía

1. C. Schabereiter-Gurtner *et al.* Novel real-time PCR assay for detection of *Helicobacter pylori* infection and simultaneous clarithromycin susceptibility testing of stool and biopsy specimens. *Journal of Clinical Microbiology* 2004; 42(10): 4512-4518.
2. J.G. Kusters *et al.* Pathogenesis of *Helicobacter pylori* infection. *Clinical Microbiology Reviews* 2006; 19(3): 449-490.
3. F. Can *et al.* Urease activity and urea gene sequencing of coccoid forms of *Helicobacter pylori* induced by different factors. *Current Microbiology* 2008; 56(2): 150-155.
4. R.X. Tang *et al.* Diversity of *Helicobacter pylori* isolates in expression of antigens and induction of antibodies. *World Journal of Gastroenterology* 2008; 14(30): 4816-4822.

5. S.K. Patel *et al.* Diagnosis of *Helicobacter pylori*: what should be the gold standard? *World Journal of Gastroenterology* 2014; 20(36): 12847-12859.
6. M. Varbanova *et al.* Chronic gastritis-an update. *Best Practice & Research: Clinical Gastroenterology* 2014; 28(6): 1031-1042.
7. Abadi A.T.B. *et al.* Resistance to clarithromycin and gastroenterologist's persistence roles in nomination for *Helicobacter pylori* as high priority pathogen by World Health Organization. *World Journal of Gastroenterology* 2017; 23(35): 6379-6384.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: March 2019

FOR INFORMATION PURPOSES ONLY

FOR INFORMATION PURPOSES ONLY

FOR INFORMATION PURPOSES ONLY

FOR INFORMATION PURPOSES ONLY

FOR INFORMATION PURPOSES ONLY

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC