# PANÓPTICO RÁPIDO CONCENTRADO 1:10

Para diagnóstico "in vitro"

## PRINCIPIO

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CLSI Guidelines and Standards, CLSI, Wayne, P.A

Gurr, E. (1965) "The rational use of dyes in Biology", p. 115. Leonard Hill, Londres.

Gurr, E. (1971) "Synthetic dyes in Biology, Medicine and Chemistry". Academic Press. London & New York.

Maree, L.; du Plessis, S.S.; Menkvels, R. and van der Horst, G. Human Reproduction, 25 (6), 1369 – 1382 (2010).

# QUICK PANOPTIC CONCENTRATED (1:10)

For in vitro diagnostic

## PRINCIPLE

The stains which make up the Quick Panoptic stain combine polychromy and the quality of classic haematology staining methods (May-Grünwald, Giemsa, Wright) with a very quick execution time of just 15 seconds. The technique is performed by immersion in the staining solutions.

As with the other Romanowsky stains, the basic dyes bind to the acidic components of cells, nucleic acids, granules in neutrophils and acidic proteins which are stained a relatively deep red-purple colour, while the acid dyes bind to haemoglobin, basic components of cell structures and eosinophil granules.

Its use allows the differential staining of blood cells. The result of the staining may be influenced by several factors, such as the fixation, the staining time and the pH value of the rinsing water. If the pH is too basic the staining will be more blue, and if the pH is too acidic, the staining will be more pink.

## DIAGNOSTIC USE

For the differential staining of peripheral blood and bone marrow smears. It can also be used for the staining of semen samples.

## REAGENTS

### Kit 3 x 100 mL

Contains:

A.- Reagent No. 1	1 x 10 mL	Ref. 99 03 50
B.- Reagent No. 2	1 x 100 mL	Ref. 99 03 52
C.- Reagent No. 3	1 x 100 mL	Ref. 99 03 54

## PREPARATION OF THE WORKING REAGENT:

Reagent No. 1: Dilute the vial to 1,000 mL with methanol.

Reagent No. 2: Dilute the vial to 1,000 mL with deionised water.

Reagent No. 3: Dilute the vial to 1,000 mL with deionised water.

Once diluted to 1,000 mL, the 3 reagents are ready for use.

## REAGENT COMPOSITION

A. Aqueous solution of hexamethyl-p-rosaniline

B. Buffered aqueous solution of xanthene.

C. Buffered aqueous solution of thiazine

## STORAGE AND STABILITY

Reagents which are stored at 15-30°C and protected from light are stable until the expiry date stated on the label. Containers must always be kept tightly closed. Once diluted they are stable at room temperature (15 - 30°C), for at least three months, if contamination and excessive evaporation is prevented.

Over time, a light precipitate may form in some reagents. This does not affect their functionality.

## MATERIAL REQUIRED (NOT SUPPLIED)

General-purpose laboratory material.

Slides for blood smears.

Manual or automated staining system.

Microscope.

## PRECAUTIONS

The products are for professional use. All the reagents must be handled with care by trained technicians. The safety data sheet should be consulted before using the reagents.

Waste disposal must be carried out according to the local regulations in force.

## SAMPLE

Air-dried blood smears. It is recommended that they are thin and homogeneous in order to obtain a better fixation of the dye without overstaining.

Staining should be carried out during the two hours following preparation in order to obtain good results. Old smears may stain irregularly.

The ideal sample is capillary blood, but if venous blood is used then EDTA should be used as an anticoagulant. The use of Heparin is not advised.

Handle the samples with care due to their potentially infectious nature.

## QUALITY CONTROL

Following the quality control practices defined by the CLSI (formerly NCCLS) is recommended. The staining method indicated is the one which is used routinely in our laboratories.

## PROCEDURE

### HAEMATOLOGY STAIN

1. Leave to air dry.
2. Set out stains No. 1, No. 2 and No. 3 in Wertheim (or similar) cuvettes.
3. Immerse the rack with the slides for 5 seconds (5 immersions lasting 1 second each) in staining solution No. 1 and drain.
4. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 2 and drain again.
5. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 3.
6. Rinse with tap water and leave to dry.
7. Examine under the microscope.

### SPERM CELLS' STAIN

8. Prepare a smear with 15 µl of fresh semen on a standard slide and leave to air dry for at least 10 minutes.
9. Set out stains No. 1, No. 2 and No. 3 in Wertheim (or similar) cuvettes.
10. Immerse the rack with the slides for 5 seconds (5 immersions lasting 1 second each) in staining solution No. 1 and drain.
11. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 2 and drain again.
12. Immerse for another 5 seconds (5 immersions lasting 1 second each) in staining solution No. 3.
13. Rinse with tap water and leave to dry.
14. Examine under the microscope.

## RESULTS

Red blood cells: Pale or deep pink.

Platelets: Pale violet or purple.

Neutrophils: Dark blue nucleus. Pink cytoplasm with red-violet granulations.

Eosinophils: Blue nucleus. Blue cytoplasm with red or red-orange granules.

Basophils: Dark blue or purple nucleus. Purple, almost black, granules.

Lymphocytes: Violet nucleus. Sky blue cytoplasm.

Monocytes: Very pale violet nucleus. Sky blue cytoplasm.

Head of sperm cell: Dark violet.

Acrosome of sperm cell: Violet, but a lighter shade of violet than the head.

Mid-piece and tail of the sperm cell: Dark violet.

Background: Pale pink.

## NOTES

The staining intensity is proportional to the staining time. The stain must be kept away from moisture in order to obtain optimal results. The staining intensity can be varied by modifying the number of immersions in stains No. 2 and No. 3, depending on whether it is wished to emphasise the pink or blue shades.

The cuvettes containing the stains, especially that of stain **No. 1**, should be kept covered, in order to prevent excessive evaporation, which could lead to colour deviations compared to normal staining.

Before the content of the cuvette is depleted, add new quantities of the staining solution on a daily basis in order to maintain an adequate level. Renew all the content from time to time. Using buffered water instead of running tap water is recommended for the rinsing processes. Running tap water has a varied pH and salt content, which, along with the high concentrations of chlorine, may alter the results of the staining and produce inconsistent results.

Each user may apply the different versions of this procedure, both manual and automated, adapting it to their standard method.

## REFERENCES

CLSI Guidelines and Standards, CLSI, Wayne, P.A

Gurr, E. (1965) "The rational use of dyes in Biology", p. 115. Leonard Hill, London.

Gurr, E. (1971) "Synthetic dyes in Biology, Medicine and Chemistry". Academic Press. London & New York.

Maree, L.; du Plessis, S.S.; Menkvels, R. and van der Horst, G. Human Reproduction, 25 (6), 1369 – 1382 (2010).

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## QUÍMICA CLÍNICA APLICADA S.A.