# **RIDASCREEN<sup>®</sup> Legionella**

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Enzyme immunoassay for the antigen detection of  
L. pneumophila serogroup 1 in urine

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---

RIDA<sup>®</sup> and RIDASCREEN<sup>®</sup>  
are registered trademarks of R-Biopharm AG  
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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## Contents

|   | page |
|---|------|
| 1. Intended use .....   | 18   |
| 2. General .....  | 18   |
| 3. Test principle .....   | 19   |
| 4. Reagents provided .....                                      | 20   |
| 5. Reagents required but not provided .....                     | 20   |
| 6. Warnings and precautions for the users .....                 | 21   |
| 7. Storage instructions .....                                   | 21   |
| 8. Indication of instability or deterioration of reagents ..... | 22   |
| 9. Specimen collection and storage .....                        | 22   |
| 10. Test procedure .....  | 23   |
| 11. Analysis .....  | 25   |
| 12. Remarks about the test procedure and interpretation .....   | 27   |
| 13. Clinical results .....                                      | 27   |

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# RIDASCREEN<sup>®</sup> Legionella

Enzyme immunoassay for the antigen detection of  
*L. pneumophila* serogroup 1 in urine

## 1. Intended use

The RIDASCREEN<sup>®</sup> Legionella kit is an in vitro diagnosticum for the detection of *Legionella pneumophila* serogroup 1 antigen in urine. It is intended as an adjunct to culture for the presumptive diagnosis of past or current Legionnaires' Disease.

## 2. General

*Legionella pneumophila* was first isolated and characterized in 1976 after a major outbreak of pneumonia at an American Legion convention in Philadelphia, Pennsylvania; thus the name, " Legionnaires' Disease". *Legionella pneumophila* is now recognized to be a common cause of community-acquired and nosocomial pneumonia. *Legionella pneumophila* serogroup 1 causes the majority of cases of legionellosis (80 – 85%). The disease presents with a broad spectrum of illness, ranging from a mild, often non-productive cough and fever to complete respiratory and multi-organ failure. In addition, gastrointestinal symptoms are prominent in Legionella pneumonia. Legionnaires' Disease can be acquired by the inhalation of aerosols associated with air handling systems, respiratory therapy equipment, and whirlpool baths. The elderly are seen as most susceptible, although children and neonates are also affected. Risk factors include cigarette smoking, immunosuppression, alcohol consumption and concomitant pulmonary disease. Approx. 80% of Legionella pneumonia patients excrete soluble Legionella antigen in their urine. This presents the opportunity for rapid detection of Legionella urinary antigen (LUA) in a urine specimen. Rapid diagnosis and early initiation of appropriate antimicrobial therapy can significantly reduce the mortality associated with *Legionella pneumophila*, which has a reported mortality rate of up to 25% in untreated immunocompetent patients.

Traditional laboratory methods for the detection of pneumonia caused by *Legionella pneumophila* infection require an adequate respiratory specimen (e.g. sputum, bronchial washing, transtracheal aspirate, lung biopsy) or paired serum specimens (acute and convalescent) for accurate diagnosis. Laboratory analysis of respiratory specimens includes traditional culture and fluorescent antibody staining. Alternatively, indirect fluorescent antibody testing of paired specimens can be performed on serum. Accurate diagnosis of *Legionella pneumophila* using these methods depends upon proper specimen collection and the level of technical expertise available. The characteristic lack of productive sputum by Legionnaires' Disease patients necessitates invasive procedures to obtain an adequate respiratory specimen for laboratory analysis. Furthermore, patient compliance in obtaining the specimens can be poor, and serological methods are retrospective in nature.

Because of these limitations, urine is the preferred specimen for collection, transport and detection in all phases of the disease. Berdal originally reported the presence of a specific soluble Legionella antigen in the urine of patients with Legionnaires' Disease. Subsequently, immunoassays were developed for the detection of this soluble Legionella antigen in urine. The RIDASCREEN® Legionella test employs the ELISA technique for the detection of soluble Legionella urinary antigen (LUA) in patients with *Legionella pneumophila* serogroup 1 infection.

### **3. Test principle**

This test is an enzyme immunoassay (ELISA). On the surface of the microtiter wells, specific antibodies against *Legionella pneumophila* serogroup 1 are bound. Urine samples and controls are pipetted into the wells. A second monoclonal antibody conjugated to horseradish peroxidase is added and then incubated at 37 °C. The simultaneous incubation results in the Legionella urinary antigen being sandwiched between the solid phase and enzyme-linked antibodies. Unbound POD-conjugate is removed by washing.

Substrate A and Substrate B are added to the wells and incubated. The enzyme bound in the wells converts the colorless Substrate to a blue color. Addition of Stop Solution converts the color from blue to yellow. The absorption is measured at 450 nm wavelength. The color intensity is directly proportional to the amount of antigen present in the sample.

Alternatively, the results can be visually read using the Visual Interpretation Card provided.

## 4. Reagents provided

The reagents in one package are sufficient for 96 determinations.

Each test kit contains:

- 1x 12 Microtiterstrips with 8 wells each (divisible) in a frame;  
coated with purified antibodies (rabbit) against  
*Legionella pneumophila* serogroup 1
- 1 x Wash Concentrate (100 ml; 20x conc.)  
contains Thimerosal
- 1 x Positive Control (1.5 ml);  
human urine containing *Legionella pneumophila* serogroup 1 antigen;  
ready to use, contains preservative
- 1 x Negative Control (1.5 ml);  
human urine;  
ready to use, contains preservative
- 1 x Conjugate (6.5 ml);  
HRP-conjugated, purified antibody (rabbit) against  
*Legionella pneumophila* serogroup 1;  
dyed red, ready to use
- 1 x Substrate A (12 ml);  
Tetramethylbenzidine (TMB); ready to use
- 1 x Substrate B (12 ml);  
Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>); ready to use
- 1 x Stop Solution (12 ml);  
1 M phosphoric acid
- 1 x Plate lid;  
for covering stripwells during incubations to prevent evaporation
- 1 x Visual Interpretation Card
- 1 x Instructions for use

## 5. Reagents required but not provided

### 5.1. Reagents

- Distilled or deionized water

## 5.2. Accessories

- micropipet for volumes of 50 µl and 100 µl
- microplate washer or multichannel pipet, alternatively wash bottle
- graduated cylinder to dilute Wash Solution (1:20)
- incubator (non-CO<sub>2</sub>!) or heat block capable to 37 °C
- paper towels
- timer with alarm

### optional:

- microplate reader (450 nm)

## 6. Warnings and precautions for the users

The Positive Control inside the kit contains autoclaved antigen of *Legionella pneumophila* serogroup 1. However, the Positive Control as well as the patient samples should be considered potentially contagious and be treated with the necessary safety precautions.

The Conjugate and the Wash Concentrate contain Thimerosal. Contact with skin or mucous membranes must be avoided.

Hydrogen peroxide can cause cauterization. Handle with care!

The Stop Solution contains 1 M phosphoric acid. Avoid contact with skin and clothing!

All reagents and materials coming in contact with potential infectious specimens must be treated with disinfectants or autoclaved at 121 °C for at least one hour.

An exchange of individual reagents between kits of different lot numbers is not possible.

Microtiterstrips and reagents must not be used if pouch is damaged or vials are leaking.

## 7. Storage instructions

All reagents have to be stored at 2 – 8 °C and can be used up to the expiry date printed on the labels. A quality warranty cannot be given beyond the kit expiration date.

Stop Solution as well as Wash Concentrate may be stored at 2 – 30 °C.

The diluted Washing Buffer has a shelf life of one week if stored at 2 – 30 °C.

Allow reagents and Microwell Strips to get room temperature before use. To avoid moisture within the strips, do not take the strips out of the foil bag before having

reached room temperature. The foil bag should be opened with a pair of scissors without detaching the fastener. Return any unused strips to the foil bag, reseal and store them directly at 2 – 8 °C.

The colorless Substrate-solutions must be protected from exposure to direct light to avoid deterioration or coloration by autoxidation. If the Substrate-solutions turn blue, the reagents should be discarded.

## **8. Indication of instability or deterioration of reagents**

The following criteria may indicate a reagent deterioration:

- a turbidity or a blue coloration of the Substrate-solutions prior to their use
- an absorbance value of the Negative Control higher than 0.120
- an absorbance value of the Positive Control lower than 4 x O.D. of Negative Control

If using visual interpretation, the Negative Control must be clearly negative (colorless), and the Positive Control should be +2 or greater, showing a definite yellow color.

Reagents which show a turbidity must not be used in the assay.

## **9. Specimen collection and storage**

Urine specimens should be collected in standard sterile containers. Whenever possible, urine specimens should be shipped in leak-proof containers at 2 – 8 °C or frozen.

Urine samples can be used fresh or frozen. Fresh samples that have not been preserved should be stored at 2 - 8 °C or at room temperature and should be tested within 24 h.

Alternatively, specimens may be stored at 2 – 8 °C for up to 14 days or frozen (-70 °C) for longer periods before testing. Deep freezing does not pose a negative influence on the test results. Repeated thawing and freezing must be avoided.

Urine specimens stored below 8°C which contain excess urates, phosphates or other dissolved salts may develop salt crystals after storage. Insure all samples are at room temperature prior to running in the assay. For dissolving the crystals the samples can be heated up to 37 °C. Contaminated or very turbid samples should be filtered before use.

## 10. Test procedure

### 10.1. Preliminary comments

The test has to be used only by experienced laboratory personnel. Please refer to guidelines for safety regulations in medical laboratories. The test protocol must be followed strictly.

Bring all reagents and the microtiterstrips to room temperature before use. Mix the reagents well before use. Reproducibility in any EIA depends on exact pipetting, the observance of incubation times and temperature and the consistency of wash sequences. During the washing steps, take care that all wells are filled with buffer and that the liquid is completely removed from the wells. Do not allow microwells to dry between steps.

In addition, do not use wash programs that include a soaking step.

Avoid direct sunlight during all incubations. Covering the microtiter plate is recommended.

Except of the Washing Buffer, all reagents are ready to use.

### 10.2. Preparation of the Washing Buffer

1 part of the Wash Concentrate is diluted with 19 parts of distilled water. The diluted Washing Buffer has a shelf life of one week if stored at 2 – 8 °C.

### 10.3. Preparation of the samples

Urine specimens are used directly in the assay. It has to be insured that all samples are at room temperature prior to running in the assay. Contaminated or very turbid samples should be filtered before use.

#### **Remark:**

**Urine specimens stored below 8°C which contain excess urates, phosphates or other dissolved salts may develop salt crystals after storage. For dissolving the crystals the samples can be heated up to 37 °C.**

#### **Attention!**

**Positive Control as well as Negative Control are ready for use and have to be tested in the same way as the patients' samples.**

#### 10.4. First incubation

After a sufficient number of wells for samples and controls have been placed into the microtiter well holder, 100 µl of the Positive Control, the Negative Control and the samples are pipetted into separate wells. 50 µl of conjugate are added to each well. After covering the microwells with the enclosed plate lid, mix by gently swirling on tabletop and incubate at 37 °C for 50 minutes.

#### 10.5. Washing

Remove plate lid and decant or aspirate all wells into a waste container with a disinfectant. Ensure complete removal of the liquid from the wells by tapping the inverted plate onto absorbent paper. Fill 400 µl of prepared Washing Buffer in all wells. Repeat the wash cycle 4 times. Be sure to remove residual washing solution by firmly tapping the inverted microwells on absorbent paper after final washing.

**If a microplate washer is used, be sure that the liquid is completely sucked off. After final washing step the inverted microwells should be firmly tapped on absorbent paper.**

#### 10.6. Second incubation

Add 50 µl of Substrate A and 50 µl of Substrate B into each well \*). Cover the plate with the enclosed plate lid and incubate the plate for 10 min at 37 °C. Following the incubation, the reaction is stopped by adding 100 µl Stop Solution to each well. After careful mixing (soft tapping on the edge of the plate) the absorbance is measured at 450 nm against an air blank or read visually by using the Visual Interpretation Card.

\*) alternatively, pre-mix immediately before use adequate equal volumes of Substrate A and Substrate B in a clean trough, mix gently and pipette 100 µl into all wells.

#### **Summary of the test procedure**

1. Bring all reagents to room temperature
2. Dilute the Wash Concentrate
3. Pipet 100 µl of the samples, the Positive and Negative Control into the microwells
4. Add 50 µl of Conjugate; 50 minutes incubation at 37 °C
5. Discard the incubate and wash 4 times with 400 µl of Washing Buffer

6. Add 50 µl of Substrate A and 50 µl of Substrate B; 10 minutes incubation at 37 °C
7. After addition of 100 µl Stop Solution visual or spectrophotometric determination (450 nm)

## 11. Analysis

### 11.1. Quality control

For quality control, the Positive and Negative Control must be included in each assay, to ensure reagent stability and correct performance of the assay procedure. Alternatively, known positive and negative urine specimens may be run as controls. Controls are intended to monitor for substantial reagent failure. The assay run is correct, if the OD for the enclosed Negative Control is below 0.120 and the OD for the enclosed Positive Control is above the four-fold value of the Negative Control. If the Negative Control yields absorbance values >0.120 this may indicate insufficient washing. If the expected control values are not fulfilled please check the following before repeating the test:

- Expiration date of the reagents
- Calibration of the used instruments
- Exact test procedure
- Visual examination of kit components for signs of contamination, deterioration or leakage; Substrate solutions must not be used if turned blue

If the control data are not fulfilled after repeating, please contact your local distributor of R-Biopharm.

Do not report patient test results if the Controls do not meet above criteria for acceptability.

### 11.2. Spectrophotometric determination

The cut-off is determined by multiplying the mean value of the Negative Control by 4.

$$\text{cut-off} = 4 \times \text{absorbance value of the Negative Control}$$

Samples are considered **positive** if the absorbance value is higher than the determined cut-off.

Samples are considered **negative** if the absorbance value is lower than the determined cut-off.

Example:

Absorbance value Negative Control = 0.052  
Absorbance value Positive Control = 1.986  
Absorbance value sample 1 = 1.253

Cut-off calculation: cut-off =  $4 \times 0.052 = 0.208$

Result sample 1:  $1.253 > 0.208 = \text{positive}$

### 11.3. Visual Determination

A sample with an at least definite yellow color (+2) is considered **positive** for Legionella urinary antigen.

A sample with no yellow color or faint yellow color (+1) is considered **negative** for Legionella urinary antigen.

Please refer to the Visual Interpretation Card included in this kit.

Negative result = - = no yellow color  
+1 = faint yellow color

Positive result = +2 = definite yellow color  
+3 = bright yellow color  
+4 = intense yellow color

### 11.4. Results reporting

| Test result                 | Recommended report   |
|-----------------------------|--|
| <b>Positive</b> (> Cut-off) | presumptive positive for the presence of <i>L. pneumophila</i> serogroup 1 antigen in urine, suggesting current or past infection.   |
| <b>Negative</b> (< Cut-off) | presumptive negative for <i>L. pneumophila</i> serogroup 1 antigen in urine, suggesting no recent or current infection. Legionnaires' disease cannot be ruled out since other serogroups and species may also cause disease. |

## 12. Remarks about the test procedure and interpretation

The RIDASCREEN® Legionella assay detects *Legionella pneumophila* serogroup 1 antigen in urine specimens. A relation between the absorbance value and the clinical relevance is not given. Assay results should always be interpreted in connection to the clinical diagnosis.

A positive result does not exclude the presence of other pathogens. A negative result does generally not exclude a Legionellosis infection.

Excretion of Legionella antigen in urine may vary depending on the individual patient and the stage of their disease. Some individuals have been shown to excrete antigen for an extended period of time, so a positive ELISA reaction may reflect a recent but not active infection.

Early treatment with appropriate antibiotics may also decrease antigen excretion in some individuals. Antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards.

## 13. Clinical results

### 13.1. Expected values

266 randomly selected, presumed negative, uncharacterized urine specimens were tested using the RIDASCREEN® Legionella ELISA. The specimens were from four reference laboratories and included both fresh and frozen samples. 263 specimens tested negative for Legionella urinary antigen. The three presumed negative specimens tested low positive. The overall accuracy was 98.9 %.

| O.D. value<br>450 nm | 0,000-0,099  | 0,100-0,149 | 0,150-0,199 | 0,200-0,249 | 0,250-0,299 | > 0,299                         |
|----------------------|--------------|-------------|-------------|-------------|-------------|---------------------------------|
| number of<br>samples | 250<br>(94%) | 7<br>(2,6%) | 4<br>(1,5%) | 1<br>(0,4%) | 1<br>(0,4%) | 3<br>(1,1%)<br>(false positive) |

### 13.2. Sensitivity and Specificity

The performance of the RIDASCREEN® Legionella ELISA was evaluated in a clinical study by a leading infectious disease reference laboratory. This was a retrospective study on a panel of clinically well-characterized urine specimens. During the study, 234 urine specimens were tested with both the RIDASCREEN® Legionella ELISA and a commercial Legionella Urinary Antigen

EIA. This trial consisted of urine specimens from 94 *Legionella pneumophila* serogroup 1 culture-positive individuals, and 140 culture-negative individuals (including non-*Legionella* pneumonia, bacteremia, miscellaneous pulmonary diseases, and urinary tract infections). The performance of the RIDASCREEN® *Legionella* ELISA and the commercial LUA ELISA against culture results are presented below:

|  | <b>Culture (+)*<br/>n=94</b> | <b>Culture (-)**<br/>n=140</b> |
|--|------------------------------|--------------------------------|
| <b>RIDASCREEN® <i>Legionella</i> ELISA for LUA (+) n=104</b> | 89                           | 15                             |
| <b>RIDASCREEN® <i>Legionella</i> ELISA for LUA (-) n=130</b> | 5                            | 125                            |
| <b>Other EIA-Test for LUA (+) n=105</b>                      | 82                           | 23                             |
| <b>Other EIA-Test for LUA (-) n=129</b>                      | 12                           | 117                            |

\* Positive respiratory culture for *Legionella pneumophila* serogroup 1 (n=94)

\*\* Negative respiratory, blood or urine culture for *Legionella pneumophila* serogroup 1 (n=140)  
total n=234

The performance characteristics for the RIDASCREEN® *Legionella* ELISA Test from this data set are as follows:

Sensitivity = 89/94 = 94.7 % (95% Confidence interval = 90.1 – 99.2%)

Specificity = 125/140 = 89.3 % (95% Confidence interval = 84.2 – 94.4%)

Accuracy = 214/234 = 91.4 % (95% Confidence interval = 84.2 – 94.4%)

A second, prospective, clinical trial correctly identified as negative 127 freshly obtained urine specimens from non-*Legionella* patients with lower respiratory infections, sepsis, or urinary tract infections and abnormal urine specimens.

|  | <b>Culture (+)*<br/>n=0</b> | <b>Culture (-)**<br/>n=127</b> |
|--|-----------------------------|--------------------------------|
| <b>RIDASCREEN® <i>Legionella</i> ELISA for LUA (+) n=0</b>   | 0                           | 0                              |
| <b>RIDASCREEN® <i>Legionella</i> ELISA for LUA (-) n=127</b> | 0                           | 127                            |

\* Positive respiratory culture for *Legionella pneumophila* serogroup 1 (n=0)

\*\* Negative respiratory, blood or urine culture for *Legionella pneumophila* serogroup 1 (n=127)  
total n=127

### 13.3. Cross reactivity

In the clinical trial, of the 140 frozen Legionella-negative specimens, 50 had pneumonia (non-Legionella), 25 were bacteremic (non-Legionella), 15 had other pulmonary conditions, and 50 had urinary tract infection. Of the 15 specimens that produced a false positive result, 14 were from patients with bacteremia (*Streptococcus pneumoniae*), and one patient had a urinary tract infection (*Escherichia coli*, *Enterobacter cloacae*). However, in a second clinical trial for specificity using freshly obtained urine specimens, none of the patients with sepsis or pneumonia (23), none of the patients with a positive urine culture (25) and none of the patients with urinary tract infection (9) produced a false positive result.

### 13.4. Precision

The Intra-Assay Variation of the RIDASCREEN® Legionella ELISA was determined on three urine specimens containing different antigen concentrations (negative, low positive, high positive). The urine specimens were run in 30 replicates on three kit lots. Each specimen was correctly identified as being positive or negative 98.5 % of the time.

| <b>sample</b> | <b>mean O.D. value</b> | <b>Standard Deviation</b> | <b>Coefficient of Variation (%)</b> |
|---------------|------------------------|---------------------------|-------------------------------------|
| negative      | 0.08                   | 0.011                     | 13.8                                |
| low positive  | 0.34                   | 0.034                     | 10.0                                |
| high positive | 1.31                   | 0.100                     | 7.6                                 |

The Inter-Assay Variation of the RIDASCREEN® Legionella ELISA was determined using 8 urine specimens containing different antigen concentrations. Specimens were made by spiking a pool of normal human urines with heat extract from *Legionella pneumophila* at concentrations that produced low positive (B, F), mid positive (A), high negative (D) and negative (G, H) O.D. results. These eight samples were masked and distributed to three sites for evaluation. Each site tested all eight spiked urine specimens in triplicate, using three different kit lots, on three different days. The specimens were correctly identified as being positive or negative in 100 % of the cases.

| <b>sample</b> | <b>mean O.D. value</b> | <b>Standard Deviation</b> | <b>Coefficient of Variation (%)</b> | <b>n =</b> |
|---------------|------------------------|---------------------------|-------------------------------------|------------|
| A             | 1.72                   | 0.23                      | 13.31                               | 26         |
| B             | 1.07                   | 0.16                      | 14.79                               | 26         |
| C             | 2.89                   | 0.27                      | 9.42                                | 21*        |
| D             | 0.15                   | 0.02                      | 16.12                               | 26         |
| E             | 2.65                   | 0.32                      | 12.19                               | 25*        |
| F             | 0.74                   | 0.13                      | 17.18                               | 26         |
| G             | 0.08                   | 0.01                      | 16.18                               | 26         |
| H             | 0.09                   | 0.02                      | 21.48                               | 26         |

\* Sample data was excluded from statistical analyses because O.D. value exceeded the threshold of plate reader instrument. Samples above threshold were considered positive for Legionella pneumophila serogroup 1 urinary antigen.

## Appendix

### Literature

1. Fraser D.W., Tsai T.R., Orenstein W., Parkin W.E., Beecham H.J., Sharrar R.G., Harris J., Mallison G. F., Martin S.M., McDade J.E., Shepherd C.C., Brachman P.S., and The Field Investigation Team. 1977. Legionnaires' disease: Description of an epidemic of pneumonia. *N. Engl. J. Med.* 29:1189-1197.
2. Stout, J.E., Yu. V.L., 1977. Legionellosis. *N. Eng. J. Med.* 337:682-687.
3. Kohler, R.B., Wilde III C., Johnson W., Joly J., Wheat L.J., Baker R., and Missett M. 1988. Immunologic diversity among serogroup 1 *Legionella pneumophila* urinary antigens demonstrated by monoclonal antibody enzyme-linked immunosorbent assays. *J. Clin. Microbiol.* 26:2059-2063.
4. Kohler, R.B. Winn W.C., Jr., and Wheat L.J. 1984. Onset and duration of urinary antigen excretion in Legionnaires' disease. *J. Clin. Microbiol.* 20:605-607.
5. Berdal, B.P., Farshy C.E., and Feeley J.C. 1979. Detection of *Legionella pneumophila* antigen in urine by enzyme-linked immunospecific assay. *J. Clin. Microbiol.* 9:575-578.
6. Tilton, R.C. 1979. Legionnaires' disease antigen detected by enzyme-linked immunosorbent assay. *Ann. Intern. Med.* 90:697-698.
7. Kohler, R.B., Zimmerman S.E., Wilson E., Allen S.D., Edelstein P.H., Wheat L.J., and White A. 1981. Rapid radioimmunoassay diagnosis of Legionnaires' Disease. *Ann. Intern. Med.* 94:601-605.
8. Bibb, W.F. P.M. Arnow, Thacker L., and McKinney R.M. 1984. Detection of soluble *Legionella pneumophila* antigens in serum and urine specimens by enzyme-linked immunosorbent assay with monoclonal and polyclonal antibodies. *J. Clin. Microbiol.* 20:478-482.
9. Tang, P.W., and Toma S.. 1986. Broad-spectrum enzyme-linked immunosorbent assay for detection of *Legionella* soluble antigens. *J. Clin. Microbiol.* 24:556-558.