## **RIDASCREEN® Leishmania Ab**

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Enzyme immunoassay for the detection of  
antibodies against Leishmania spec.

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RIDA® and RIDASCREEN®  
 are registered trademarks of R-Biopharm AG  
 Manufacturer: R-Biopharm AG, Darmstadt, Germany

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## Contents

	page
1. Intended use.....	13
2. General.....	13
3. Test principle.....	14
4. Reagents provided.....	15
5. Materials required but not provided.....	15
6. Warnings and precautions for the users.....	16
7. Storage instructions.....	16
8. Indication of instability or deterioration of reagents.....	17
9. Specimen collection and storage.....	17
10. Test procedure.....	17
11. Analysis.....	19
12. Remarks about the test procedure and interpretation.....	20

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## RIDASCREEN® Leishmania Ab

Enzyme immunoassay for the detection of antibodies against Leishmania spec.

### 1. Intended use

The RIDASCREEN® Leishmania Ab test is an enzyme immunoassay (EIA) for the detection of specific antibodies against Leishmania (species: donovani, major, mexicana, braziliensis) in human serum.

### 2. General

Leishmaniasis is a parasitic disease caused by different species of a flagellate protozoa genus Leishmania. The most characteristic clinical forms of the disease are:

- Visceral leishmaniasis or kala-azar: produced by *L. donovani* (India, Africa), *L. infantum* (Mediterranean zone, Asia, Middle East) and *L. chagasi* (South America). The phagocytosed parasites are numerous in the reticuloendothelial (RE) cells of the spleen, liver, lymph nodes, bone marrow, intestinal mucosa and other organs. These organs tend to enlarge due to the enormous growth of the RE cells. Incubation periods last 1 to 4 months. The initial fever is intermittent and rise to 40 °C. A marked leucopenia with relative monocytosis and lymphocytosis, anemia and thrombocytopenia develop, rendering the patient especially susceptible to secondary infections. Untreated disease usually progresses to a fatal termination within 2 years, although fulminating infections may cause death within few weeks.
- Cutaneous leishmaniasis: Produced by *L. major* or *tropica* (Old World) and by *L. mexicana* and *braziliensis* (America). The disease is known with different names in each country. In humans the disease is limited to the cutaneous tissues. Incubation periods last from 2 weeks to several months. It is the least severe form of the disease.

- c) American or mucocutaneous leishmaniasis: Mainly produced by *L. braziliensis*, although other species produced similar mucocutaneous lesions. The clinical appearance and histopathology of american leishmaniasis are identical to cutaneous leishmaniasis, except that it may produce later mucous membrane involvement. The mucosal lesions are painful and can cause great deformity with erosion of the nasal septum, palate, or larynx.

The life cycle involves an alternate existence in a vertebrate and an insect host. The natural host, besides humans, are rodents and carnivores. The invertebrate hosts are sandflies of the genus *Phlebotomus* and *Lutzomyia* (America). After parasites within infected cells are ingested with the blood meal, they transform within the insects from the amastigote (without flagella) to the promastigote form (with flagella) and multiply in the gut. In 8 to 20 days, the anterior gut and pharynx are partially blocked by flagellates. When the sandfly attempts a subsequent blood meal, some of the infective promastigotes are dislodged and introduced into the skin. Transmission can also occur by contamination of wounds.

Leishmaniasis has a world wide distribution, mainly in tropic and subtropic zones. The incidence of the disease is rising and infection is calculated to threaten 350 million people. It is one of the six most important infectious diseases in the world.

### 3. Test principle

On the surface of the microtiter wells, inactivated soluble antigens of *Leishmania* are bound. Diluted serum samples, controls and calibrator are pipetted into the wells and incubated at room temperature. Present antibodies bind to the immobilized antigens. Unbound material is removed in a washing step.

In a second step, an anti-human Ig-conjugate (peroxidase-conjugated IgG antibody) is added, which binds to the antibodies bound before. After incubation, unbound conjugate is removed by washing. Substrate (urea peroxide) and chromogen (TMB) are added to the wells and incubated at room temperature. The enzyme bound in the wells converts the colorless substrate/chromogen to a blue color. Addition of stop solution converts the color from blue to yellow. The absorption is measured at 450 nm wavelength (optional: reference wavelength  $\geq 620$  nm).

## 4. Reagents provided

### 4.1. Test kit for 96 determinations (Art. No. K 7121)

The reagents in one package are sufficient for 96 determinations. Each test kit contains:

- 1 x 12 Microwell Strips (breakable) with 8 wells each in a frame; coated with inactivated *Leishmania* antigen; in a resealable foil bag
- 1 x Sample Buffer (50 ml, pH 7.4), ready to use, dyed yellow, colorless lid
- 1 x Washing Buffer (50 ml, pH 7.4), 20x conc., dyed blue, brown lid
- 1 x Calibrator (1.2 ml), green screw cap; inactivated human serum, ready to use
- 1 x Positive Control (1.2 ml), red screw cap; inactivated human serum, ready to use
- 1 x Negative Control (1.2 ml), colorless screw cap; inactivated human serum, ready to use
- 1 x Anti-human Ig-Conjugate (12 ml), red cap; HRP-conjugated IgG-antibody (goat), ready to use
- 1 x Substrate (6 ml), green cap; urea peroxide, ready to use
- 1 x Chromogen (6 ml), blue cap; tetramethylbenzidine (TMB), ready to use
- 1 x Stop Solution (6 ml), yellow cap; 1 M sulfuric acid
- 1 x Instructions for use

### 4.2. Test kit for 48 determinations (Art. No. K 7122)

The reagents in one package are sufficient for 48 determinations. Each test kit contains:

- 1 x 6 Microwell Strips (breakable) with 8 wells each in a frame; coated with inactivated *Leishmania* antigen; in a resealable foil bag

Other kit components are identical to Art. No. K 7121 (see point 4.1.).

## 5. Reagents required but not provided

### 5.1 Reagents

- Distilled or deionized water

## 5.2 Accessories

- Test tubes
- Vortex mixer
- Micropipets for volumes of 10 - 100 µl and 100 - 1000 µl
- Measuring cylinder (1000 ml)
- Microplate washer or multichannel pipet
- Microplate reader (450 nm, optional: reference wavelength  $\geq$  620 nm)
- Absorbent paper

## 6. Warnings and precautions for the users

The control sera and the calibrator have been tested for HIV- and HCV-Ab as well as for HBsAg and were found to be negative. However, they as well as the patient samples should be considered potentially contagious and be treated with the necessary safety precautions.

Hydrogen peroxide can cause cauterization. Handle with care!

The stop solution contains 1 M sulfuric acid. Avoid contact with skin and clothing!

All reagents and materials coming in contact with potential infectious specimens must be treated with disinfectants or autoclaved at 121 °C for at least one hour.

An exchange of reagents between kits of different lot numbers is not possible.

## 7. Storage instructions

All reagents have to be stored at 2 - 8 °C and can be used up to the expiry date printed on the labels. Microbial contamination has to be avoided. A quality warranty cannot be given beyond the kit expiration date.

The diluted washing buffer has a shelf life of 4 weeks if stored at 2 - 8 °C.

Allow reagents and microwell strips to get room temperature before use. To avoid moisture within the strips, do not take the strips out of the foil bag before having reached room temperature. The foil bag should be opened with a pair of scissors without detaching the fastener. Return any unused strips to the foil bag, reseal and store them directly at 2 - 8 °C.

The colorless chromogen must be protected from exposure to direct light to avoid deterioration or coloration by autoxidation. If the chromogen has turned blue, the reagent should be discarded.

## 8. Indication of instability or deterioration of reagents

The following criteria may indicate a reagent deterioration:

- a turbidity or a blue coloration of the colorless chromogen prior to its use
- an absorbance value (O.D.) of the negative control at 450 nm  $>$  0.3
- an absorbance value (O.D.) of the positive control at 450 nm  $<$  0.6

## 9. Specimen collection and storage

The RIDASCREEN® Leishmania Ab EIA has been evaluated for the investigation of human serum samples. Repeated freezing and thawing of the samples as well as microbial contamination must be avoided. The application of heat treated, lipemic, hemolytic, icteric or turbid samples can lead to wrong results.

The sample material can be stored for up to 1 week at 2 - 8 °C if the test cannot be carried out immediately. A prolonged storage of samples is possible at -20 °C.

If they have been stored at 2 - 8 °C, diluted samples can be used for up to 7 hours.

## 10. Test procedure

### 10.1. Preliminary comments

Bring all reagents and the microwell strips for at least 30 minutes to room temperature before use. Mix the reagents well before use. Reproducibility in any EIA depends on exact pipetting, the observance of incubation times and temperature and the consistency of wash sequences.

During the washing steps, take care that all wells are filled with buffer and that the liquid is completely removed from the wells. Do not allow microwells to dry between steps.

Avoid direct sunlight during all incubations. Covering the microtiter plate is recommended.

Except the washing buffer, all reagents are ready to use.

### 10.2. Preparation of the washing buffer

1 part of the concentrated washing buffer is diluted with 19 parts of distilled water. Crystals in the buffer concentrate can be dissolved in a waterbath at 37 °C.

Add 50 ml of the concentrated washing buffer to a 1000 ml graduated cylinder. Bring the final volume to 1000 ml with distilled or deionized water. The diluted washing buffer has a shelf life of 4 weeks if stored at 2 - 8 °C.

### 10.3. Preparation of the samples

Before starting the test, serum samples have to be diluted 1:20 with the sample buffer.

e.g. 10 µl serum + 190 µl sample buffer

#### **Attention!**

**The controls included in the kit are ready to use and must not be diluted.**

### 10.4. First incubation

After insertion of a sufficient number of cavities into the microwell holder, 100 µl of the positive control, the calibrator, the negative control and the diluted sera are pipetted into the wells and incubated for 10 min at room temperature. It is recommended to use the calibrator in double determination in each run.

### 10.5. Washing

Decant or aspirate all microwells into a waste container with a disinfectant. Ensure complete removal of the liquid from the microwells by tapping the inverted plate onto absorbent paper. Then wash all wells 5 times with 300 µl of prepared washing buffer. Be sure to remove residual washing solution by firmly tapping the inverted microwells on absorbent paper after single washing steps.

**If a microplate washer is used, take care that the washer is adjusted to the used microplate type.**

### 10.6. Second incubation

Add 100 µl or 2 drops of anti-human Ig-conjugate to all wells. Incubate the plate for 5 min at room temperature.

### 10.7. Washing

Wash 5 times according to step 10.5.

### 10.8. Third incubation

Add 50 µl or 1 drop of substrate and 50 µl or 1 drop of chromogen into each well. Incubate the plate for 10 min at room temperature in the dark. Following the incubation, the reaction is stopped by adding 50 µl or 1 drop of stop solution to each well. After careful mixing (soft tapping on the edge of the plate) the absorbance is measured in a microplate reader at 450 nm (optional: reference wavelength  $\geq$  620 nm).

#### **Remark:**

**Highly positive patient samples can cause dark precipitates of the chromogen.**

#### **Summary of the test procedure**

1. Bring all reagents to room temperature
2. Dilute the washing buffer
3. Dilute the serum samples
4. Pipet 100 µl of the positive control, the calibrator (double determination), the negative control and the diluted samples into the microwells; 10 minutes incubation at room temperature
5. Discard the incubate and wash 5 times with 300 µl of washing buffer
6. Add 100 µl or 2 drops of anti-human Ig- conjugate; 5 minutes incubation at room temperature
7. Discard the incubate and wash 5 times with 300 µl of washing buffer
8. Add 50 µl or 1 drop each of substrate and chromogen; 10 min incubation at room temperature in the dark
9. After addition of 50 µl or 1 drop stop solution spectrophotometric determination at 450 nm (optional: reference wavelength  $\geq$  620 nm)

## **11. Analysis**

### 11.1. Quality control

For the quality control, positive control, calibrator (double determination) and negative control must be carried along with each test run. The test was carried out correctly, if the positive control shows an absorbance value (O.D.) at 450 nm greater than 0.6. The ratio between the O.D. of the positive control and the O.D. of the calibrator must be greater than 1.5. The absorbance value of the negative control at 450 nm must show a value lower than 0.3. Additionally, the O.D. of the negative control must be less than the O.D. of the calibrator.

### 11.2. Calculation of the sample ratio

- a) Calculate the average O.D. of the calibrator. This value is the cut-off value of the assay.
- b) Divide the sample O.D. by the value obtained in point a).

e.g.: calibrator well 1            O.D. = 0.541  
         calibrator well 2        O.D. = 0.554  
         sample                    O.D. = 1.441

cut-off value = Fehler!= 0.548

sample ratio = Fehler!= 2.63

### 11.3. Test result

Tab. 1: Valuation of sample ratio

	negative	equivocal	low positive	middle positive	high positive
sample ratio	< 0.9	0.9 - 1.1	1.1 - 1.5	1.5 - 2.5	> 2.5

### 12. Remarks about the test procedure and interpretation

The RIDASCREEN® Leishmania Ab EIA detects specifically antibodies against Leishmania (species: donovani, major, mexicana, braziliensis). It should be used in suspected cases of Leishmaniasis. Assay results should always be interpreted in connection to the clinical diagnosis and other diagnostic results.

Equivocal and low positive results can happen in these cases:

- Persons that were bitten by infected Phlebotomus at endemic areas and did not develop the disease but have a low antibody level.
- Persons that presently have or suffered from asymptomatic disease, so they still have a low antibody level.
- Persons that suffered from the disease in the past and still have a low antibody level.
- Persons that are developing the disease.
- Unidentified and unspecific reaction factors.

These cases should be further investigated, retesting a fresh new sample and including clinical symptoms and additional tests. In case the same result is obtained, the test must be repeated with a new sample after 2-4 weeks.

Middle and high positive results can happen in persons that are suffering from the disease or have suffered a recent infection, or with past infections that still have a high antibody level. It is necessary to compare these results with clinical symptoms and additional tests.

Negative antibody findings cannot exclude a leishmaniasis. Due to low antibody titer at early time of an infection, the test can show negative results. If a clinical suspicion subsists, after two to four weeks another patient's sample should be tested.

It has been founded that in South America, this test may cross react in some occasions with serum positive to T.cruzi (or Chagas disease). Therefore, special caution must be taken in areas of South America where Chagas disease is present, when using this test.

## Appendix

### Literature

1. Bray, R.S. and Lainson, R., Trans. Roy. Soc. Trop. Med. Hyg., 60, 605-609 (1966).
2. Bray, R.S. and Lainson, R., Trans. Roy. Soc. Trop. Med. Hyg., 61, 490-505 (1967).
3. Bray, R.S., Ecol. Dis., 4, 257-267 (1982).
4. Hommel, M., Peters, W., Ranque, J., Quilici, M. and Lanotte, G., Ann. Trop. Med. Parasitol., 72, 213-218 (1978).
5. Kager, P.A., Rees, P.H., Welde, B.T., Hockmeyer, W.T. and Lyerly, W.H., Trans. R. Soc. Trop. Med. Hyg., 75, 556-559 (1981).
6. Lainson, R., Trans. Roy. Soc. Trop. Med. Hyg., 77, 569-596 (1983).
7. Low-A-Chee, R.M., Rose, P. and Ridley, D.S., Ann. Trop. Med. Parasitol., 77, 255-260 (1986).
8. Ranque, J. and Quilici, M., Journal of Parasitology, 56, 277-278 (1970).