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



Parvovirus

recomLine
Parvovirus B19
IgG [Avidity], IgM

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Suspicious of infection < 6 – 8 weeks				Sum of points						Test result	
				IE1	CM2	p150	p85	gB1	gB2		
				G	M	G	M	G	M	G	M
				1	1	5	3	6	4	1	2
				1	1	1	1	1	1	1	1
IgG	React Control	IgG	IgM	IE1	CM2	p150	p85	gB1	gB2	7 points	positive
Avidity										≥2 avidity antigens with low avidity.	suspect. inf. < 14 weeks
IgM										8 points	positive
Suspicious of infection < 14 weeks											
IgG										14 points	positive
Avidity										≥2 avidity antigens with low avidity.	suspect. inf. < 14 weeks
IgM										4 points	positive
Suspicious of infection > 12 weeks											
IgG										14 points	positive
Avidity										p150 and ≥ 1 avidity antigen with high avidity.	suspect. inf. > 12 weeks
IgM										0 points	negative
Suspicious of infection > 24 weeks											
IgG										7 points	positive
Avidity										-	no additional statement
IgM										0 points	negative

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day after symptoms	recomWell Parvovirus B19			recomLine Parvovirus B19							interpretation Line	total interpretation (incl. information of Line)	
	IgG	IgM	interpretation ELISA	Cutoff	VP-2p	VP-N	VP-1S	VP-2r	VP-C	NS-1			result
2	12 U/ml negative	221 U/ml positive	suspected acute infection	IgG	+/-	+/-	+/-	+/-	-		negative	suspected acute infection	< 2 weeks (IgG seroconversion) avidity: no relevance, because IgG negative
				avidity							low avidity*		
				IgM	2+	+	+/-	+	+		positive		
8	35 U/ml positive	225 U/ml positive	suspected acute infection	IgG	+	+	+	2+	+		positive	suspected acute infection	< 26 weeks (VP-C-IgG positive, while at the same time IgM is positive)
				avidity							low avidity		
				IgM	2+	2+	+	2+	+		positive		
15	186 U/ml positive	177 U/ml positive	suspected acute infection	IgG	3+	3+	3+	3+	3+		positive	suspected acute infection	< 26 weeks (VP-C-IgG positive, while at the same time IgM is positive)
				avidity							intermediate avidity		
				IgM	2+	2+	+	2+	+		positive		
43	242 U/ml positive	70 U/ml positive	suspected acute infection	IgG	3+	3+	3+	3+	2+	+/-	positive	status after infection (weeks to months)	< 26 weeks (VP-C-IgG positive, while at the same time IgM is positive) avidity: > 4 weeks
				avidity							high avidity		
				IgM	+	+	+/-	2+	+		positive		
71	305 U/ml positive	60 U/ml positive	suspected acute infection	IgG	2+	3+	3+	3+	2+	+	positive	status after infection (weeks to months)	> 6 weeks and < 26 weeks (NS-1-IgG positive and VP-C-IgG positive, while at the same time IgM is positive) avidity: > 4 weeks
				avidity							high avidity		
				IgM	+	+	-	+/-	+/-		positive		

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