NanoEnTek

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※ Product code : 17712 (30Tests/Kit)

▶ INTENDED PURPOSE

BioTracer™ TB Rapid Card is intended for the detection of antibodies to *Mycobacterium tuberculosis* in human serum or plasma.

▶ EXPLANATION OF THE TEST

Tuberculosis is a highly infectious disease caused by *Mycobacterium tuberculosis*. The methodologies routinely used for diagnosis of tuberculosis are chest X-ray examination and culture technique. Chest X-ray cannot early detect the tuberculosis and culture technique require a minimum of 2 weeks before a definitive diagnosis of tuberculosis. Serology test can overcome these disadvantages with its simple and rapid method. **BioTracer™ TB Rapid Card** is a chromatographic immunoassay kit for rapid qualitative determination for *M.tuberculosis* infection using human serum or plasma specimen.

▶ PRINCIPLE OF THE METHOD

BioTracer™ TB Rapid Card is based on the principle of an immuno-chromatography in vitro test for the qualitative determination of antibodies against *M.tuberculosis* in human serum or plasma. When the specimens are added to the sample pad, they move to the conjugate pad and resuspend the *M.tuberculosis* antigen conjugated gold complex that is dried on the conjugate pad. The mixtures move along the membrane by capillary action and react with the *M.tuberculosis* antigen that are immobilized on the test reaction zone. If antibodies against *M.tuberculosis* are present enough in the sample, a colored band in the test reaction zone will be appeared. If there are no antibodies against *M.tuberculosis* or not sufficient in the sample, the area will remain colorless. The mixture continues to move to the control reaction zone and forms a red or purple color, indicating the test is working properly and the result is valid.

▶ CONTENTS

1. Test device
2. Dropper
3. Package insert

▶ SPECIMEN COLLECTION AND PREPARATION

1. Serum or plasma specimens should be tested immediately after collection.
2. Do not leave the specimens at room temperature for prolonged period. Specimens may be stored at 2-8°C for up to 2 days. For long term storage, specimens should be kept below -20°C.

※ Specimens containing precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

▶ TEST PROCEDURE

Allow the specimens and sealed pouch containing the test devices to room temperature prior to use.

1. Remove the test device from the sealed pouch, and place it on a clean and flat surface.
2. Apply 3 drops (about 120 µL) of serum or plasma specimens into the sample well (S) of the test device.
3. Read the test result at 10~15 minutes.

※ Caution: Perform the test immediately after removing the test device from the foil pouch.

▶ READING AND INTERPRETATION OF RESULT

1. Control (C) band means that the test is working properly.
2. Test (T) band indicates the test result.

NEGATIVE:

The presence of only Control(C) band indicates a negative result.

POSITIVE:

The presence of Test band (T) with Control (C) band indicates a positive result.

INVALID:

If control (C) band is not appeared in the result window after performing test, the result is considered invalid.

※ The directions may not have been followed correctly or the test device may have been deteriorated. It is recommended that the specimens be re-tested with the new device.

▶ FOLLOW-UP ACTION

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

▶ LIMITATIONS OF THE METHOD

A negative result does not preclude the possibility of infection with *M. tuberculosis*. Other clinically available tests are required if questionable results are obtained.

▶ PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

The sensitivity for *M.tuberculosis* was more than 82.5% (66/80). The specificity was 92.7% (306/330).

TB Patients	BioTracer™ TB Rapid Card		Total
	Positive	Negative	
Positive	66	14	80
Negative	24	306	330
Total	90	320	410

2. Reproducibility

- 1) Within run performance test was determined by one analyst with ten times for 3 different positive and 1 negative control specimens. There were no variation within the test devices.
- 2) Between run performance test was determined by three analyst with 3 different lots for 3 different positive and 1 negative control specimens. There were no variation between different analyst.

3. Interference test

No interference was found with bilirubin (10mg/dL), hemoglobin (20mg/dL) or triglycerides (600mg/dL) on the sensitivity and specificity of the **BioTracer™ TB Rapid Card**.

▶ WARNING AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Do not eat or smoke while handling specimens.
3. Wear protective gloves while handling specimens. Wash hands thoroughly afterward.
4. Do not use test kit if the packing is damaged or the seal is broken.
5. Avoid splashing or aerosol formation while handling specimens.
6. Clean up spilled specimens thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, tested kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.

▶ STORAGE AND SHELF-LIFE

BioTracer™ TB Rapid Card should be stored at 2~30°C (36~86°F). The test device is sensitive to humidity as well as to heat. Do not use it beyond the expiration date, 24 months from manufacturing date.

▶ LITERATURE REFERENCES

1. Ian A.C, Oumou B.S. Pulmonary tuberculosis: diagnosis and treatment. BMJ. 2006; 332: 1194-1197.
2. Manju R.P, Tehmina M., Harald G.W. 2007. Immunohistochemical diagnosis of abdominal and lymph node tuberculosis by detecting *Mycobacterium tuberculosis* complex specific antigen MPT64. Diagnostic Pathology. 2007; 2: 36.
3. Liselotte V.A, Miranda L., Robert Z., Raymond P.B. Tuberculosis risk varies with the duration of HIV infection: a prospective study of European drug users with known date of HIV seroconversion. AIDS. 2003; 17(8): 1201-1208.

Manufactured by
NanoEnTek Inc.
 851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,
 Gyeonggi-do, 18531, Korea
 Tel : +82-2-6220-7942 / Fax : +82-2-6220-7999
 ivdst@nanoentek.com
 www.nanoentek.com

MT Promed Consulting GmbH
 Altenhofstrasse 80, 66386
 St. Ingbert, Germany