VIASURE

Real Time PCR Detection Kits

by **CerTest**
BIOTEC

Tropical Panel I

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VIASURE *Tropical Panel I* Real Time PCR Detection Kit 12 x 8-well strips, low profile

VS-TP0112L

VIASURE *Tropical Panel I* Real Time PCR Detection Kit 12 x 8-well strips, high profile

VS-TP0112H

VIASURE *Tropical Panel I* Real Time PCR Detection Kit 96-well plate, low profile

VS-TP0113L

VIASURE *Tropical Panel I* Real Time PCR Detection Kit 96-well plate, high profile

VS-TP0113H

ENGLISH

1. Intended use

VIASURE *Tropical Panel I* Real Time PCR Detection Kit is designed for the specific and differentiation of Zika (ZIKV), Dengue (DENV), Chikungunya (CHIKV), West Nile (WNV), Yellow fever (YFV) and/or Mayaro (MAYV) viruses in clinical samples from patients with signs and symptoms of Zika, Dengue, Chikungunya, West Nile, Yellow fever and/or Mayaro viruses infection. This test is intended for use as an aid in the diagnosis of the ZIKV, DENV, CHIKV, WNV, YFV and MAYV in combination with clinical and epidemiological risk factors. RNA is extracted from specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for Zika, Dengue, Chikungunya, West Nile, Yellow fever, and Mayaro viruses.

2. Summary and Explanation

Zika (ZIKV), Dengue (DENV) and Chikungunya (CHIKV) are arthropod-borne viruses (arboviruses) that share mosquitos of the *Aedes* genus as common vector, specifically *A. aegypti* and *A. albopictus*. Environmental changes, urbanization and the globalization of travel have enhanced the spread of ZIKV, DENV and CHIKV resulting in epidemics and the co-circulation in overlapping geographic areas as well as the possibility of viral co-infection within a single host. Currently the major endemic regions include Africa, Southeast Asia, the Pacific Islands and the Americas.

Besides, these arboviruses cause similar clinical presentations, especially in the initial stages of infection (fever, headache, skin rash, muscle and joint pain, and gastrointestinal symptoms). In fact, due to neither virus possesses any specific distinguishing clinical features and the outcomes and management strategies for these three viruses are vastly different, an early and accurate diagnosis is imperative.

Therefore, ZIKV, DENV and CHIKV must be initially included in the differential diagnosis for a patient with suspicious clinical symptoms who is living in or returning from travel to an endemic area. Laboratory diagnostic tests are based on the detection of the virus, viral components (antigens or nucleic acid), or the host immunologic response to the virus. However, due to the cross-reactivity of the antibodies of these arboviruses limits the use of serology, Real Time RT-PCR is a detection method commonly used during the acute phase of the infection. In order to do that, clinical specimens that can be tested including blood, serum, plasma, urine and others.

West Nile virus (WNV) is considered the most important causative zoonotic pathogen of viral encephalitis worldwide. WNV belongs to the genus *Flavivirus* and is a member of the Japanese encephalitis serogroup. Phylogenetic analysis has revealed two main lineages of WNV, lineage 1, which is widely distributed and further subdivided in three clades (A, B and C) and lineage 2. While several additional lineages have been proposed during the last years.

In nature, there is an enzootic transmission cycle of WNV between birds (hosts) and mosquitoes (bridge vectors) of the genus *Culex* (mainly *Cx. pipiens*, *Cx. perexiguus* and *Cx. modestus*) and *Aedes*. But this virus can infect other vertebrates and humans, between which it has been reported additional methods of transmission: blood transfusion, organ transplantation, transplacental transmission, and via breastmilk. Most human infections with WNV (~80%) are asymptomatic or characterized by a mild-febrile syndrome and flu-like malaise. Nevertheless,

fewer than 1% infection may progress to more severe neuroinvasive disease (WNND) and even death, depending on the host immune response and the viral strain. WNND usually encompasses three different syndromes: encephalitis, meningitis, and acute flaccid paralysis.

Laboratory diagnosis is generally accomplished by testing of serum or cerebrospinal fluid (CSF) to detect WNV-specific IgM antibodies and confirmed by plaque-reduction neutralization tests (PRNTs). However, the cross-reactivity with infections caused by other flaviviruses remains the major problem with most serological diagnostic tests. In this regard, different PCR-based protocols have been developed to detect viral RNA and can be performed on serum and CSF, both collected early in the course of illness, and additional tissue specimens as urine, where it can be detected much longer and at higher concentrations.

Yellow fever virus (YFV) is an RNA virus that belongs to the *Flaviviridae* family. The "yellow" in the name refers to the jaundice caused by liver involvement that affects some patients. Yellow fever is an acute viral haemorrhagic zoonotic disease transmitted between humans and from monkeys to humans by mosquitoes, primarily by *Aedes (aegypti)* and *Haemogogus*.

The incubation period is generally 3–6 days from infection until illness. Many people do not experience symptoms, but when these do occur during acute febrile phase patients could develop myalgia, headache, back pain, anorexia, nausea, and sometimes vomiting; symptoms typically resolve within 1 week. Small proportions of patients enter a second, more toxic phase within 24 hours of recovering from initial symptoms, and develop severe symptoms (with high fever, jaundice, bleeding, and kidney damage). Half of the patients who enter this toxic phase die; due to there is currently no specific anti-viral drug. Although, specific care to treat dehydration, liver and kidney failure, and fever improves outcomes. Therefore, vaccination and mosquito control are the most important means of preventing yellow fever.

Yellow fever is difficult to diagnose, especially during the early stages. It can be confused with severe malaria, leptospirosis, viral hepatitis (especially the fulminating forms of hepatitis B and D), other hemorrhagic fevers, infection with other flaviviruses (Dengue, Zika and West Nile viruses), and poisoning.

Usually, the detection of yellow fever virus is carried out by virus cultivation on susceptible cells followed by immunofluorescence, ELISA, immunohistochemical examination or PCR. For quantifying the load of infectious particles in cell culture or serum samples, plaque assay is still the commonly used method. However, this is cumbersome and time-consuming technique. With the recent development of Real-Time PCR, a rapid and accurate alternative to quantify viral load in blood is available.

Mayaro virus (MAYV) is an enveloped, single-stranded positive RNA virus which belongs to the genus *Alphavirus*, family *Togaviridae*. Its transmission cycle is thought to occur mainly through mosquito vectors, especially those of genus *Haemagogus*, but *Aedes* spp. mosquitoes may also be competent vectors. These vectors transmit the virus to primates and birds and sometimes to humans.

Phylogenetic studies using a fragment of the structural polyprotein open reading frame suggest the classification of MAYV in 2 different genotypes: D (includes isolates from all countries) and L (contains strains detected only in Brazil).

The disease caused by Mayaro lasts for 3-5 days and includes fever, headache, rash, myalgia, arthralgia and, occasionally, arthritis, which can be very debilitating and can persist for months. This symptomatology is quite similar to other arthropod-borne virus, being sometimes clinically indistinguishable from Dengue, Chikungunya fever, malaria, rabies, measles or other arboviral diseases, so it is required to use a reliable and specific test to avoid this misdiagnosis.

This Real Time PCR Detection Kit enable a reliable diagnosis of Mayaro virus by the amplification of the *nsp1* gene, which has been shown to be the best target for detection due to its highly conserved sequences. Because MAYV is blood-borne and mosquito-transmitted, blood is the usual human sample collected. Serum and plasma are prepared from this and used to determine Mayaro infection.

3. Principle of the procedure

VIASURE *Tropical Panel I* Real Time PCR Detection Kit is designed for the diagnosis of the Zika, Dengue, Chikungunya, West Nile, Yellow fever and/or Mayaro viruses in clinical samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of specific genes using specific primers and a fluorescent-labelled probes (Table 1, Pathogens detected in each reaction tube and target genes).

Tropical Panel I 8-well strip contains in each reaction well the following reaction mixes for the detection of the specific targets (Table 1):

	Code	Controls	Reaction mix placed into each well Pathogens and target genes
1	ZDC	IC	Zika, Dengue & Chikungunya
			Zika (<i>envelope gene</i>), Dengue (3'Non coding region) and Chikungunya (<i>NSP1 gene</i>)
2	WNV	IC	West Nile Virus
			West Nile Virus (a conserved sequence of genomic region 5'UTR)
3	FEV	IC	Yellow Fever
			Yellow Fever (a conserved sequence of genomic region 5'UTR)
4	MYV	IC	Mayaro Virus
			Mayaro Virus (<i>nsp1 gene</i>)

Table 1. *Tropical panel I* 8-well strips provided in VIASURE *Tropical panel I* Real Time PCR Detection Kit. Reaction mix placed into each well, pathogens detected and target genes. IC: Internal control. Note that the first reaction mix is placed into the ends of the strip.

VIASURE *Tropical Panel I* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Tropical Panel I* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to monitor PCR inhibition. Each RNA targets are amplified and detected in specific channels (FAM, HEX, ROX, and/or Cy5) and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2)(Table 5).

4. Reagents provided

VIASURE *Tropical Panel I* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 2 and Table 3:

Reagent/Material	Description	Colour	Amount
<i>Tropical Panel I</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Tropical Panel I</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Tropical Panel I* Real Time PCR Detection Kit with Ref. VS-TP0112L and VS-TP0112H.

Reagent/Material	Description	Color	Amount
<i>Tropical Panel I</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Tropical Panel I</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 3. Reagents and materials provided in VIASURE *Tropical Panel I* Real Time PCR Detection Kit with Ref. VS-TP0113L and VS-TP0113H.

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Tropical Panel I* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE *Tropical Panel I* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-TP0113L and VS-TP0113H). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.

- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- **Note that the first reaction mix is placed into the ends of the strip.**
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.

8. Test procedure

8.1. RNA extraction

Perform the sample preparation following the manual system protocol or according to the recommendations appearing in the instructions for use of extraction kit used.

For RNA extraction from clinical samples (blood, serum, plasma, urine, tissues and others) you can use your manually or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- QIAamp Viral RNA Mini Kit (Qiagen).
- NucleoMag® Pathogen (Macherey Nagel).
- NucleoSpin® RNA Virus (Macherey Nagel).
- Invisorb® Spin Universal Kit (Stratec).
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- NX-48 Viral NA Kit, using the Nextractor® NX-48 system (Genolution).
- MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.).
- ZP02003 MagPurix Viral Nucleic Acid Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.)
- ZP02013 MagPurix Viral RNA Extraction kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.)

8.2. Lyophilized positive control

Tropical Panel I Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized

Tropical Panel I Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required strips, including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plate or strips. **Note that the first reaction mix is placed into the ends of the strip (Table 1) and each strip includes two *Tropical Panel I* reactions, wells 1 to 4: ZDC, WNV, FEV and MYV reaction mixes).**

- 1) Reconstitute the number of strips you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted *Tropical Panel I* Positive Control (red vial) or Negative Control (violet vial) in the 8 wells of each strip and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler.

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM, ROX, Cy5 and/or HEX (JOE or VIC) channels following Table 5. Depending on the equipment used select the proper detection channel (see Annex 2).

The reaction mixes placed in each well of *Tropical Panel I* strip allow the detection of the specific target pathogens in the following channels (Table 5).

	Code	Channels			
		FAM	HEX	ROX	Cy5
1	ZDC	Dengue Virus	Internal Control	Chikungunya virus	Zika virus
2	WNV	West Nile virus	Internal Control	-	-
3	FEV	Yellow Fever	Internal Control	-	-
4	MYV	Mayaro Virus	Internal Control	-	-
4	MYV	Mayaro Virus	Internal Control	-	-
3	FEV	Yellow Fever	Internal Control	-	-
2	WNV	West Nile virus	Internal Control	-	-
1	ZDC	Dengue Virus	Internal Control	Chikungunya virus	Zika virus

Table 5. Pathogens identified in each detection channel.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative controls of each well from *Tropical Panel I* reaction(4 wells) and the presence of signal for *Tropical Panel I* positive control of each well from *Tropical Panel I* reaction (4 wells). Check Internal Control signal (ZDC, WNV, FEV and MYV reaction mixes) to verify the correct functioning of the amplification mix.. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer’s instructions.

A sample is considered positive if the Ct value obtained is less than 40 and the Internal Control shows or not an amplification signal. Sometimes, the detection of internal is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.


A sample is considered negative, if the sample shows no amplification signal in the detection system but the Internal Control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control (if available).

Using the following table read and analyze the results:

	Reaction Mix	Pathogens	Channels			
			FAM	HEX	ROX	Cy5
1	ZDC	Dengue Virus	Positive			
		Chikungunya virus			Positive	
		Internal Control		Positive/Negative		
		Zika virus				Positive
2	WNV	West Nile virus	Positive			
		Internal Control		Positive/Negative		
3	FEV	Yellow Fever	Positive			
		Internal Control		Positive/Negative		
4	MYV	Mayaro Virus	Positive			
		Internal Control		Positive/Negative		

Table 6. Sample interpretation. Positive: Amplification curve. Empty: No amplification curve.


Table 7. Correct run of negative and positive controls run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative controls or absence of signal in the positive wells. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells, we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with RNA extracted from serum and urine samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Zika, Dengue, Chikungunya, West Nile, Yellow fever and/or Mayaro viruses either by samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Tropical Panel I* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

VIASURE *Tropical Panel I* Real Time PCR Detection Kit (for **Zika, Dengue & Chikungunya Reaction mix**) was evaluated with samples from INSTAND and QCMD programs. The clinical performance of VIASURE assay was tested using 102 clinical specimens. The results were compared with the final EQA Reports:

	EQA Reports			
		+	-	Total
Zika, Dengue & Chikungunya Reaction mix	+	25	0	25
	-	1	76	77
	Total	26	76	102

Table 8. Comparative results for Zika Virus

	EQA Reports			
		+	-	Total
Zika, Dengue & Chikungunya Reaction mix	+	27	0	27
	-	3	72	75
	Total	30	72	102

Table 9. Comparative results for Dengue Virus

	EQA Reports			
		+	-	Total
Zika, Dengue & Chikungunya Reaction mix	+	29	1	30
	-	2	70	72
	Total	31	71	102

Table 10. Comparative results for Chikungunya Virus

VIASURE *Tropical Panel I* Real Time PCR Detection Kit (for **West Nile Virus**) was evaluated with QCMD 2014, 2016, 2017 and 2018 panels from West Nile Virus RNA EQA Programme and with INSTAND e.V.'s Virus Genome Detection – West Nile Virus 2016 and 2017. These panels consist of 64 clinical specimens dissolved in transport medium. The results were compared with those presented by EQA programme final reports. All West Nile Virus (47/64) samples could be detected. Samples contained different dilutions of West Nile Virus strains NY99 (lineage 1), Ug37 and Heja (lineage 2). In addition, WNV Negative and Non-WNV flaviviruses samples (which include Japanese Encephalitis virus/Dengue virus types 1-4/ Tick borne Encephalitis virus/ Yellow Fever virus/Zika virus) could be confirmed as a negative.

The clinical performance of VIASURE *Tropical Panel I* Real Time PCR Detection Kit (for **Yellow fever virus**) was evaluated using Tropical Virus panels (Chikungunya, Dengue, Zika and West Nile Viruses) from QCMD for several years (2014, 2015, 2016 and 2017). These QCMD panels included 9 positive samples for yellow fever virus. The results were compared with the programs final reports. Other Tropical Virus including Dengue Virus Type 1, Dengue Virus Type 2, Dengue Virus Type 3, Dengue Virus Type 4, Zika Virus, Chikungunya Virus, Japanese Encephalitis virus, West Nile Virus (strains NY99 and Heja), and Tick borne Encephalitis virus, were not detected.

The results show a high sensitivity and specificity to detect Zika, Dengue, Chikungunya, West Nile and Yellow fever viruses using VIASURE *Tropical Panel I* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Tropical Panel I* Real Time PCR Detection Kit has a detection limit of ≥ 10 RNA copies per reaction for Zika, Dengue, Chikungunya, West Nile, Yellow fever and/or Mayaro.

12.3. Analytical specificity

The specificity of the Zika, Dengue, Chikungunya, West Nile, Yellow fever and/or Mayaro viruses assays were confirmed by testing a panel consisting of different microorganisms representing the most common arbovirus. No cross-reactivity was detected between any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing			
Zika virus strain MR 766 (Uganda)	-/+	Chikungunya virus strain S27 Petersfield	-/+
Zika Virus strain 11474/16 (French Polynesia)	-/+	St Louis Encephalitis virus strain 17D	-
Zika Virus International Standard (French Polynesian, 11468/16)	-/+	West Nile virus strain H160/99	-/+
Zika Virus (African)	-/+	West Nile virus Heja	-/+
Zika virus strain PF13/251013-18 (Asian)	-/+	West Nile virus Ug37	-/+
Dengue 1 virus strain Hawaii	-/+	Yellow Fever virus strain 17D	-/+
Dengue 2 virus strain New Guinea C	-/+	Japanese Encephalitis virus	-
Dengue 3 virus strain H87	-/+	Plasmodium falciparum	-
Dengue 4 virus strain H241	-/+	Trypanosoma cruzi	-

Table 11. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Tropical Panel I* Real Time PCR Detection Kit for **Zika, Dengue & Chikungunya Reaction mix**, was evaluated against Zika virus (Zika virus MR 766 (Uganda, 1947), Zika Virus strain 11474/16 (French Polynesia), Zika Virus strain 11468/16(French Polynesia), Zika Virus (African) and Zika virus strain PF13/251013-18 (Asian), Dengue Virus (Dengue 1 virus strain Hawaii, Dengue 2 virus strain New Guinea C, Dengue 3 virus strain H87 and Dengue 4 virus strain H241) and Chikungunya Virus (Chikungunya virus S27 Petersfield and Chikungunya virus Martinique isolate (Asian genotype)), showing positive results.

The reactivity of VIASURE **West Nile Virus Reaction mix** was evaluated against West Nile virus H160/99, Heja,Ug37 and NY99, showing positive results.

The reactivity of VIASURE **Yellow Fever Virus Reaction mix** was evaluated against Yellow Fever virus strain 17D, showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁴⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁴⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁴⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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

13. Bibliography/Bibliografía

1. R.S. Lanciotti *et al.* Genetic and serologic properties of Zika virus associated with an epidemic, Yap State, Micronesia, 2007. *Emerging Infectious Diseases* 2008; 14(8): 1232-1239.
2. O. Faye *et al.* Quantitative real-time PCR detection of Zika virus and evaluation with field-caught mosquitoes. *Virology Journal* 2013; 10: 311.
3. A.C. Gourinat *et al.* Detection of Zika virus in urine. *Emerging Infectious Diseases* 2015; 21(1): 84-86.
4. D. Musso *et al.* Detection of Zika virus in saliva. *Journal of Clinical Virology* 2015; 68: 53-55.
5. D. Mussi *et al.* Potential Sexual Transmission of Zika Virus. *Emerging Infectious Diseases* 2015 21(2): 359-361.
6. J. Mlakar. Zika Virus Associated with Microcephaly. *New England Journal of Medicine* 2016.
7. C. Drosten *et al.* Rapid detection and quantification of RNA of Ebola and Marburg viruses, Lassa virus, Crimean-Congo hemorrhagic fever virus, Rift Valley fever virus, dengue virus, and yellow fever virus by real-time reverse transcription-PCR. *Journal of Clinical Microbiology* 2002; 40(7): 2323-2330.
8. A. Dumoulin *et al.* Pan-dengue virus detection by PCR for travelers returning from the tropics. *Journal of Clinical Microbiology* 2008; 46(9): 3104-3106.
9. J. Liu *et al.* Development of a TaqMan Array Card for Acute-Febrile-Illness Outbreak Investigation and Surveillance of Emerging Pathogens, Including Ebola Virus. *Journal of Clinical Microbiology* 2016; 54(1): 49-58.
10. S.K. Mardekian and A.L. Roberts. Diagnostic Options and Challenges for Dengue and Chikungunya Viruses. *BioMed Research International* 2015: 834371.
11. K.A. Tsetsarkin *et al.* Interspecies transmission and chikungunya virus emergence. *Current Opinion in Virology* 2016; 16: 143-150.
12. R.S. Lanciotti *et al.* Chikungunya virus in US travelers returning from India, 2006. *Emerging Infectious Diseases journal* 2007; 13(5): 764-767.
13. Centers for Disease Control and Prevention. (<http://www.cdc.gov/>).
14. World Health Organization. (<http://www.who.int/>).
15. C. Chancey *et al.* The global ecology and epidemiology of West Nile virus. *BioMed Research International* 2015; 2015:376230.
16. A. Rizzoli *et al.* The challenge of West Nile virus in Europe: knowledge gaps and research priorities. *Eurosurveillance* 2015; 20(20): 21135.
17. S. Linke *et al.* Detection of West Nile virus lineages 1 and 2 by real-time PCR. *Journal of Virological Methods* 2007; 146(1-2): 355-358.
18. Centers for Disease Control and Prevention. West Nile Virus (<https://www.cdc.gov/westnile/>).
19. World Health Organization. Zika virus and potential complications (<http://www.who.int/mediacentre/factsheets/fs354/en/>).
20. J. Liu *et al.* Development of a TaqMan Array Card for Acute-Febrile-Illness Outbreak Investigation and Surveillance of Emerging Pathogens, Including Ebola Virus. *Journal of Clinical Microbiology* 2016; 54(1): 49-58.
21. A.G. Fernandes-Monteiro *et al.* New approaches for the standardization and validation of a real-time qPCR assay using TaqMan probes for quantification of yellow fever virus on clinical samples with high quality parameters. *Human Vaccines & Immunotherapeutics* 2015; 11(7): 1865-1871.
22. World Health Organization. Yellow fever. <http://www.who.int/mediacentre/factsheets/fs100/en/>

23. World Health Organization. Emergencies preparedness, response. Yellow fever. <http://www.who.int/csr/disease/yellowfev/en>
24. Centers for Disease Control and Prevention. Yellow fever. <https://www.cdc.gov/yellowfever/>.
25. M.T. Oliveira Mota *et al.* Complete Genome Sequence of Mayaro Virus Imported from the Amazon Basin to Sao Paulo State, Brazil. *Genome Announcements* 2015; 3(6): e01341-15.
26. J. Lednicky *et al.* Mayaro Virus in Child with Acute Febrile Illness, Haiti, 2015. *Emerging infectious diseases* 2016; 22(11): 2000-2002.
27. M. Llagone-Barets *et al.* A case of Mayaro virus infection imported from French Guiana. *Journal of Clinical Virology* 2016; 77: 66-68.
28. M. L. Garcia de Figueiredo *et al.* Emerging alphaviruses in the Americas: Chikungunya and Mayaro. *Revista de Sociedade Brasileira de Medicina Tropical* 2014; 47(6): 677-683.
29. B. Friedrich-Jänicke *et al.* Genome analysis of Mayaro virus imported to Germany from French Guiana. *Emerging infectious diseases* 2014; 20(7): 1255-1257.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: November 2019

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC