

MLS Portable Laboratory Autoclaves

105°C to 126°C | 45°C to 135°C

Convenient sterilisation on demand offers great flexibility

(MLS-3020U)

(MLS-3751L, MLS-3781L)



MLS-3020U 47 L



MLS-3751L 50 L



MLS-3781L 75 L

Model No.	MLS-3020U	MLS-3751L	MLS-3781L
Temperature Control Range	105°C to 126°C	45°C to 135°C	
Effective Internal Capacity	47 litres	50 litres	75 litres
External Dimensions (W x D x H)*	440 x 550 x 1050 mm	478 x 632 x 748 mm	478 x 632 x 965 mm
Chamber Dimensions (dia. x depth)	300 x 670 mm	370 x 415 mm	370 x 630 mm
Net Weight (Approx.)	67 kg	61 kg	71 kg

* External dimensions of main cabinet only - see dimension drawings on flyers showing handles and other external projections.

U: CE marked

Models MLS-3751L/3781L are for laboratory use.

