

# Quality Mixer



## Quality Mixer

QA Mixer is a new fully automatic mixer/tubestripper manufactured to simplify and standardize the requirements of quality control according to directive 2004/33/EG, article 6 appendix V.2.

QA Mixer works against ergonomic injuries in hand wrists and shoulders which happens during manual operation. QA Mixer strips the tube and mix the bag automatically. QA Mixer is adjustable from 3-10 stripping cycles and 3-10 mixing cycles. For example the tube will be stripped 5 times and the bag will be mixed 10 times between each strip. QA Mixer is suitable for all blood bags on the market.

## Technical data:

Mains unit	inlet 24 V
Mixing cycle	11 cycles per minute
Dimensions	290 (L) x 253 (W) x 150 (H) mm.
Weight	5 kg
Max Weight	500 g
Power supply unit	inlet 230 VAC + - 10%, outlet 24 V
Symbol	CE marked

**LJUNGBERG & KÖGEL AB**