



IMMUNOASSAY Unit

STEROID HORMONES

Among the most complete panel available on the market

*since 1997, the Diametra steroid hormone menu
represents the status of the art*

- *Well proven quality*
 - *Reagents ready to use*
 - *Shelf life over 15 months*
 - *Each product is CE marked.*
- *Same procedure for the whole panel and same common reagents:
easy application on the most popular microplate auto-analyzers.*

Diametra is part of IDS group

Serum Steroid hormone test kits

- Free Testosterone is among the most used in the world (more than 1.000.000 tests sold in the 2014)

CODE	ITEM	SENSITIVITY	INCUBATION TIME	ANALYTICAL RANGE
DKO001	Cortisol	2,44 ng/mL	60+15 min	10 – 500 ng/mL
DKO002	Testosterone	0,10 ng/mL	60+15 min	0,2 – 16 ng/mL
DKO003	Estradiol	8,7 pg/mL	120+30 min	20 – 2000 pg/mL
DKO004	17 OH-Progesterone	0,05 ng/mL	60+15 min	0,2 – 20 ng/mL
DKO005	DHEA-S	0,03 µg/mL	60+15 min	0,1 – 10 µg/mL
DKO006	Progesterone	0,05 ng/mL	60+15 min	0,2 – 40 ng/mL
DKO007	Free Estriol	0,03 ng/mL	60+15 min	0,15 – 40 ng/mL
DKO008	Androstenedione	0,01 ng/mL	60+15 min	0,1 – 10 ng/mL
DKO015	Free Testosterone *	0,04 pg/mL	60+15 min	0,2 – 100 pg/mL
DKO019	Total Estriol	0,22 ng/mL	60+15 min	2 – 200 ng/mL
DKO053	Aldosterone	7,0 pg/mL	60+20 min	20 – 2000 pg/mL
DKO124	DHEA	0,1 ng/ml	60+15 min	0,5 – 30 ng/ml

Urine Steroid hormone test Kits

CODE	ITEM	SENSITIVITY	INCUBATION TIME	ANALYTICAL RANGE
DKO018	Urinary Cortisol	2,9 ng/mL	60+15 min	10 – 500 ng/mL

Saliva Steroid hormone test Kits

- The measurement of free (active) fraction of steroid means more diagnostic significance compared to serum analytes
- Non invasive sampling, it can be done by patient at home
- Easy sampling in children
- Unparalleled sensitivity
- Very useful for monitoring the hormone profile (repetitive sampling during the days)
- Application in veterinary field, sport medicine and assessment of stress monitoring (i.e. air force pilots)

CODE	ITEM	SENSITIVITY	INCUBATION TIME	ANALYTICAL RANGE
DKO020	Cortisol Saliva	0,12 ng/mL	60+15 min	0,5 – 100 ng/mL
DKO021	Testosterone Saliva	3,0 pg/mL	60+15 min	10– 1000 pg/m
DKO022	Estradiol Saliva	0,5 pg/mL	120+15 min	1 – 100 pg/mL
DKO024	DHEA-S Saliva	0,05 ng/mL	15+15 min	0,2 – 12 ng/mL
DKO025	Progesterone Saliva	3,8 pg/mL	60+15 min	10 – 5000 pg/mL
DKO026	Estriol Saliva	< 2,0 pg/mL	60+15 min	2,5 – 4000 pg/mL
DKO027	Androstenedione Saliva	5,0 pg/mL	60+15 min	20 – 1000 pg/mL
DKO120	Salivette™ (Sarsted)		100 tests	
DKO178	Estriol Saliva HS	1,48 pg/mL	ON + 15 min	2,0 – 315 pg/mL

Diametra Srl - Head Office

20090 Segrate - Milano - Italy
Via Calabria, 15
Ph +39.02.21.39.184 - Fax +39.02.21.33.354
V.A.T. - IT 12585340156
info@diametra.com

R&D and Manufacturing

06034 Spello - Perugia - Italy
Via Pozzuolo, 14
Ph +39 0742 24 851 - RFax +39 0742 316197

www.diametra.com